

The Biotechnology Education Company ®

EDVO-Kit #

371

PCR-based DNA Fingerprinting

Storage: See Page 3 for specific storage instructions

Experiment Objective:

The objective of this experiment is to perform PCR-based DNA fingerprinting on actual DNA samples that are cloned in plasmids and to understand the concept of this technology as applied to forensic science.

EDVOTEK, Inc. • 1-800-EDVOTEK • www.edvotek.com



Table of Contents

| | Page |
|---|-------------|
| Experiment Components Experiment Requirements Background Information | 3 3 5 |
| Experiment Procedures | |
| Experiment Overview and General Instructions | 8 |
| Module I: PCR Amplification of Crime Scene & Suspect DNA Module II: Separation of PCR Reactions by Electrophoresis | 10 |
| Agarose Gel Preparation | 12 |
| Conducting Electrophoresis | 16 |
| Staining and Visualization of DNA | 17 |
| InstaStain Ethidium Bromide | 18 |
| One-Step Staining and Destaining | |
| with InstaStain® Blue | 20 |
| InstaStain® Blue Cards | 21 |
| Study Questions | 23 |
| Instructor's Guidelines | 25 |
| Notes to the Instructor | 26 |
| Pre-Lab Preparations | 28 |
| Quantity Prep for Agarose Gel Electrophoresis | 29 |
| Experiment Results and Analysis | 30 |
| Study Questions and Answers | 31 |
| Appendices | |
| PCR Experimental Success Guidelines | 32 |
| PCR Using Three Waterbaths | 34 |
| Preparation and Handing of PCR Samples with Wax | 35 |
| Material Safety Data Sheets | 36 |

All components are intended for educational research only. They are not to be used for diagnostic or drug purposes, nor administered to or consumed by humans or animals.

THIS EXPERIMENT DOES NOT CONTAIN HUMAN DNA. None of the experiment components are derived from human sources.



EDVOTEK, The Biotechnology Education Company, and InstaStain are registered trademarks of EDVOTEK, Inc.. Ready-to-Load and UltraSpec-Agarose are trademarks of EDVOTEK, Inc.

EDVOTEK - The Biotechnology Education Company® • 1-800-EDVOTEK • www.edvotek.com

Component Quantities:

Experiment # 371 contains reagents to perform five sets of Polymerase Chain Reaction (PCR)* reactions (25 reactions total).

Sample volumes are very small. For liquid samples, it is important to quick spin the tube contents in a microcentrifuge to obtain sufficient volume for pipeting. Spin samples for 10-20 seconds at maximum speed.

*If you do not have a thermal cycler, PCR experiments can be conducted, with proper care, using three waterbaths. However, a thermal cycler assures a significantly higher rate of success.

Experiment Components

CONTENTS

- A Tubes with PCR reaction pellets[™] Each PCR reaction pellet[™] contains:
 - dNTP Mixture •
 - Tag DNA Polymerase Buffer •
 - Tag DNA Polymerase
 - MgCl2

| В | Primer Mix | -20°C Freezer |
|---|-----------------|---------------|
| С | 200 bp ladder | -20°C Freezer |
| D | DNA Template #1 | -20°C Freezer |
| Е | DNA Template #2 | -20°C Freezer |
| F | DNA Template #3 | -20°C Freezer |
| G | DNA Template #4 | -20°C Freezer |

REAGENTS & SUPPLIES:

- UltraSpec-Agarose[™]
- Electrophoresis Buffer (50x) •
- 10x Gel Loading Solution •
- InstaStain® Ethidium Bromide
- InstaStain® Blue
- 100 ml graduated cylinder (packaging for samples) ٠
- Microcentrifuge Tubes (0.5 ml) •
- PCR tubes (0.2 ml for thermal cyclers with 0.2 ml template)
- Wax beads (for waterbath option or thermal cyclers without heated lid)

REQUIREMENTS:

- Thermal Cycler (EDVOTEK catalog #541 is recommended) •
- Alternative Option: Three waterbaths (94°C, 45°C, and 72°C) ٠
- Horizontal gel electrophoresis apparatus •
- D.C. power supply ٠
- Balance ٠
- ٠ Microcentrifuge
- UV Transilluminator or UV Photodocumentation system
- UV safety goggles
- Automatic micropipets (5-50 µl) with tips
- Microwave, hot plate or burner
- Pipet pumps or bulbs
- 250 ml flasks or beakers
- Hot gloves
- Disposable vinyl or latex laboratory gloves
- Ice buckets and ice
- Distilled or deionized water



EDVOTEK - The Biotechnology Education Company® 1-800-EDVOTEK • www.edvotek.com FAX: (301) 340-0582 • email: edvotek@aol.com

3

Room Temperature





Deoxyribonucleic acid (DNA), present in the nucleus of every living cell, is the genetic material that acts as a blueprint for all of the proteins synthesized by that cell. In mammals, however, a large fraction of the total DNA does not encode protein and serves no obvious function. Polymor-

phic DNA refers to chromosomal regions that vary widely from individual to individual. By examining several of these regions within the genomic DNA obtained from an individual, one may determine a "DNA Fingerprint" for that individual. DNA polymorphisms are now widely used for determining paternity/maternity, kinship, identification of human remains, and the genetic basis of various diseases. The most widely used and far-reaching application, however, has been to the field of criminal forensics. DNA from both crime victims and offenders can now be definitively matched to crime scenes, often affecting the outcome of criminal and civil trials.

The beginning of DNA fingerprinting occurred in the United Kingdom in 1984, following the pioneering work of Dr. Alex Jeffreys at the University of Leicester. Analysis by Jeffreys led to the apprehension of a murderer in the first DNA fingerprinting case in September 1987. The first U.S. conviction occurred on November 6, 1987 in Orlando, FL. Since then, DNA analysis has been used in thousands of convictions. Additionally, over 100 convicted prison inmates have been exonerated from their crimes, including several death row inmates.

In 1990, the Federal Bureau of Investigation (FBI) established the **Co**mbined **D**NA Index System (CODIS), a system which allows comparison of crime scene DNA to DNA profiles in a convicted offender and a forensic (crime scene) index. A match of crime scene DNA



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved

Suspect #2 matches

Crime Scene

Background Information

EDVO-Kit #

37

PCR-based DNA Fingerprinting



to a profile in the convicted offender index indicates a suspect for the crime, whereas a match of crime scene DNA to the forensic index (a different crime scene) indicates a serial offender. CODIS has now been used to solve dozens of cases where authorities had not been able to identify a suspect for the crime under investigation.

The first step in forensic DNA fingerprinting is the collection of human tissue from the crime scene or victim. These tissues include blood, hair, skin, and body fluids. The sample, often present as a stain, is treated with a detergent to rupture (lyse) cell membranes and obtain DNA for further analysis (Figure 1). One early method, called Restriction Fragment Length Polymorphism (RFLP) analysis, involves digesting DNA with restriction enzymes, separating the fragments by agarose gel electrophoresis, transferring the DNA to a membrane, and hybridizing the membrane with probes to polymorphic regions. This method is statistically very accurate, but requires relatively large amounts of DNA and takes several days to complete. Because of the time and DNA requirements, the RFLP method is no longer used in forensics, but remains in use in certain medical genetics-based tests.

More recently, the polymerase chain reaction (PCR) has been used in forensics to analyze DNA (Figure 2). This technique requires much less (500-fold) DNA than RFLP analysis and is much less time-consuming. PCR amplification (Figure 2) uses an enzyme known as Tag polymerase. This enzyme, originally purified from a bacterium that inhabits hot springs, is stable at very high (near boiling) temperatures. Also included in the PCR reaction mixture are two (15-30 nucleotide) synthetic oligonucleotides, known as "primers" and the extracted DNA, known as the "template". The region of DNA to be amplified is known as the "target". In the first step of the PCR reaction, the template complimentary DNA strands are separated (denatured) from each other at 94°C, while the Tag polymerase remains stable. In the second



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

step, known as annealing, the sample is cooled to an intermediate temperature, usually 40°-65°C, to allow hybridization of the two primers, one to each of the two strands of the template DNA. In the third step, known as extension, the temperature is raised to 72°C and the *Taq* polymerase adds nucleotides to the primers to complete the synthesis of the new complementary strands. These three steps - denaturation, annealing, and extension - constitute one PCR "cycle". This process is typically repeated for 20-40 cycles, amplifying the target sequence exponentially (Figure 2). PCR is performed in a thermal cycler, an instrument that is programmed to rapidly heat, cool and maintain samples at designated temperatures for varying amounts of time.

In forensics, PCR is used to amplify and examine highly variable (polymorphic) DNA regions. These are regions that vary in length from individual to individual and fall into two categories: 1) variable number of tandem repeats (VNTR) and 2) STR (short tandem repeats). A VNTR is a region that is variably composed of a 15-70 base pair sequence, typically repeated 5-100 times. An STR is similar to a VNTR except that the repeated unit is only 2-4 nucleotides in length. By examining several different VNTRs or STRs from the same individual, investigators obtain a unique DNA profile for that individual which is unlike that of any other person (except for an identical twin).

In this experiment, students and teachers are encouraged to design their own crime scene scenario and come up with a plan to test their crime-solving skills by using tools such as PCR.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM EDVO-Kit #

| Experiment Overvie | ew ar | nd General Instructions | |
|---------------------------|-------------------------|--|--|
| | BE | FORE YOU START THE EXPERIMEN | IT: |
| | 1. | Read all instructions before starting the | experiment. |
| | 2. | If you will be conducting PCR using a the lid, also read the Appendix entitled "Pre Samples with Wax". | ermal cycler without a heated eparation and Handling PCR |
| | | If you will be using three waterbaths to pendices entitled "Polymerase Chain Rea and "Handling Samples with Wax Overla | conduct PCR, read the two ap- action Using Three Waterbaths" ays". |
| | 3. | Write a hypothesis that reflects the expe tal outcomes. | riment and predict experimen- |
| | EX | PERIMENT OBJECTIVE: | |
| | Th on coi | e objective of this experiment is to perform actual DNA samples that are cloned in pla ncept of this technology as applied to fore | m PCR-based DNA fingerprinting asmids and to understand the ensic science. |
| | BF | RIEF DESCRIPTION OF THE EXPERIN | IENT: |
| | In exe Stu ter | this experiment, students will conduct a Pe ercise on DNA from a simulated crime scer udents will be encouraged to constitute th PCR, students will analyze the amplified I | CR-based DNA fingerprinting ne and four different suspects. ne profiles of the individuals. Af- DNA segments on agarose gels. |
| | Th | is experiment has two modules: | |
| | I. II. | PCR Amplification of Crime Scene ar Separation of PCR Reactions by Elec | nd Suspect DNA trophoresis |
| | GI | EL SPECIFICATIONS: | |
| | Th | is experiment requires a gel with the follo | wing specifications: |
| | | Recommended Gel Size: Number of Samples Wells: Placement of the Well-former Template: Gel Concentration Required: | 7 x 14 cm (long tray) 6 First set of notches 1.0% |
| | 1 | | |



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

EDVO-Kit

Experiment Overview and General Instructions

LABORATORY SAFETY

- 1. Gloves and goggles should be worn routinely as good laboratory practice.
- 2. Exercise extreme caution when working with equipment that is used in conjunction with the heating and/or melting of reagents.
- 3. DO NOT MOUTH PIPET REAGENTS USE PIPET PUMPS.
- 4. Exercise caution when using any electrical equipment in the laboratory.
 - Although electrical current from the power source is automatically disrupted when the cover is removed from the apparatus, first turn off the power, then unplug the power source before disconnecting the leads and removing the cover.
 - Turn off power and unplug the equipment when not in use.
- 5. EDVOTEK injection-molded electrophoresis units do not have glued junctions that can develop potential leaks. However, in the unlikely event that a leak develops in any electrophoresis apparatus you are using, IMMEDIATELY SHUT OFF POWER. Do not use the apparatus.
- 6. Always wash hands thoroughly with soap and water after handling reagents or biological materials in the laboratory.





Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

Module I: PCR Amplification of Crime Scene & Suspect DNA

1. Each student group should obtain the following items from the instructor:

The PCR reaction pelletTM contains *Taq* DNA polymerase, the four deoxytriphosphates, Mg+2 and buffer. Sample volumes are very small. For liquid samples, it is important to quick spin the tube contents in a microcentrifuge to obtain sufficient volume for pipeting. Spin samples for 10-20 seconds at maximum speed.

EDVO-Kit #

37

If your thermal cycler is equipped with a heated lid, proceed directly to polymerase chain reaction cycling. If your thermal cycler does not have a heated lid, add one wax bead to the tube before proceeding to polymerase chain reaction cycling. The wax bead will melt to form a layer of oil that will protect the PCR incubation reaction from evaporation. Crime scene DNA Suspect #1 DNA Suspect #2 DNA Suspect #3 DNA Suspect #4 DNA 5 PCR beads (in tubes) Primer Mix

Label tubes containing PCR beads with the different DNAs from the list above. Also include your group designation. Label the tubes accordingly:

Crime scene PCR Suspect #1 PCR Suspect #2 PCR Suspect #3 PCR Suspect #4 PCR

SETTING UP PCR REACTIONS:

Each individual PCR reaction should be prepared as follows:

Tap the PCR reaction tube to assure that the PCR reaction pelletTM is at the bottom of the tube. Add the following to the PCR tube:

- 20 µl of Primer Mix
- 5 µl of corresponding DNA
- Gently mix the reaction tube.
- Place each sample on ice until the remaining samples are prepared.
- Note use a clean pipet tip for each individual



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

EDVO-Kit #

Module I: PCR Amplification of Crime Scene & Suspect DNA

POLYMERASE CHAIN REACTION CYCLING

Program the thermal cycler for a total of 35 cycles. Each cycle will be:

Initial Denaturation 94°C for 3 min. 35 cycles @ 94°C for 30 sec. 45°C for 30 sec. 72°C for 30 sec. Final Extension 72°C for 3 min.

After the final extension is complete, add 5 μ l of 10x Gel Loading solution to each PCR sample. Store samples on ice until ready for electrophoresis.

NOTE: If your thermal cycler has the capability, you can to link to a program to hold samples at 4°C overnight after completing the final extension.

STOP

OPTIONAL STOPPING POINT

The samples can be held in the thermal cycler at 4° C or frozen at -20° C after addition of 5 µl of 10x Gel Loading Solution until ready for electrophoresis.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM The

Experiment

EDVO-Kit #

PCR-based DNA Fingerprinting

Module II: Separation of PCR Reactions by Electrophoresis

AGAROSE GEL REQUIREMENTS FOR THIS EXPERIMENT

• Recommended gel size:

7 x 14 cm

1.0%

To achieve better resolution of the PCR products, 7 x 14 cm gels, which can be shared by several students or groups, are recommended.

• Placement of well-former template:

First set of notches

Agarose gel concentration:

Agarose Gel Preparation

PREPARING THE GEL BED

- Close off the open ends of a clean and dry gel bed (casting tray) by using rubber dams or tape.
 - A. Using Rubber dams:
 - Place a rubber dam on each end of the bed. Make sure the rubber dam fits firmly in contact with the sides and bottom of the bed.
 - B. Taping with labeling or masking tape:
 - With 3/4 inch wide tape, extend the tape over the sides and bottom edge of the bed.
 - Fold the extended edges of the tape back onto the sides and bottom. Press contact points firmly to form a good seal.
- 2. Place a well-former template (comb) in the first set of notches at the end of the bed. Make sure the comb sits firmly and evenly across the bed.







Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

The Experiment

ерvо-кіt # **371**

Module II: Separation of PCR Reactions by Electrophoresis

CASTING THE AGAROSE GEL(S)

3. Use a 250 ml flask or beaker to prepare the gel solution.

IMPORTANT

Table

Check with your instructor regarding the concentration of the buffer you are using to prepare your gel. Use the appropriate table (A.I or A. 2) below.

If preparing the gel with concentrated (50x) buffer, use Table A.1.

| F | A. I | Indivi | dual 1.0 | % | UltraSpec- | Ag | garose™ | M | Gel |
|---|-------------|-------------------|--------------------------|---|--------------------------------------|----|----------------------------|---|-------------------------|
| | Siz | ze of Gel (cm) | Amt of Agarose (g) | + | Concentrated Buffer (50x) (ml) | + | Distilled Water (ml) | = | Total Volume (ml) |
| | 7 | ′ × 7 | 0.25 | | 0.5 | | 24.5 | | 25 |
| | 7 | x 14 | 0.5 | | 1.0 | | 49.0 | | 50 |

If preparing the gel with diluted (1x) buffer, use Table A.2.

Diluted buffer is one volume of concentrated buffer to every 49 volumes of distilled or deionized water. See Table B.

| Table A.2 | Ult | Individual traSpec-Agaro | l.0% ose™ Gel |
|---------------|-------------|-----------------------------|--------------------------------|
| Size o (cn | f Gel n) | Amt of Agarose + (g) | Diluted Buffer (1x) (ml) |
| 7 x | 7 | 0.25 | 25 |
| 7 x | 14 | 0.5 | 50 |

4. Swirl the mixture to disperse clumps of agarose powder.

5. With a marking pen, indicate the level of the solution volume on the outside of the flask.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM The Experiment

| | 6. | Hea sho | at th ould | he i d ap | mixture pear clea | to disso ar (like v | lve the aga water) witl | arose por nout any | wder. The final solution v undissolved particles. |
|-----------------------|------|---|--|---------------------------------------|--|---|--|------------------------------------|--|
| | | A. | Mie | licro | ' wave m | ethod: | | | , i |
| | | | • | C | Cover the | e flask w | /ith plastic | wrap to | minimize evaporation. |
| At high altitudes, it | | | • | Н | leat the | mixture | on High f | or 1 min | ute. |
| ven to reach boiling | | | • | S a | wirl the Ill the ag | mixture arose is | and heat completel | on High y dissolv | in bursts of 25 seconds un ved. |
| mperatures. | | В. | Ho | ot p | late met | thod: | | | |
| | | | • | C ti | Cover the ion. | e flask w | vith alumin | um foil | to prevent excess evapora- |
| | | | • | H ir | leat the ng. Boil | mixture until all | to boiling the agaro | over a l se is con | burner with occasional swir npletely dissolved. |
| | | | | C n | heck the ot comple | solution etely dis | carefully. If solved. | you see " | 'crystal" particles, the agarose |
| | 7. | Coc care tior has the ma | ol th eful n of s occ e solu rkec | he a I sw f he curr lutio | agarose s virling to eat. If de red, add on up to n the fla | solution promo etectable distilled the ori sk in ste | to 60°C w te even dis e evaporat d water to ginal volur ep 6. | ith sipa- ion bring ne | DO NOT POUR BOILING HOT AGAROSE INTO THE GEL BED. |
| | Aft | er t | he g | gel | l is cool | ed to 6 | 0°C: | | Hot agarose solution may irreversibly warp the bed. |
| | lf y | ou a | are i | usi | ing rubb | oer dan | ns, go to s | tep 9. | |
| | lf y | ou a | are (| usi | ing tape | e, contir | nue with s | tep 8. | |
| | 8. | Sea tior | al the n fro | ne ir om | nterface leaking. | of the g | gel bed and | d tape to | o prevent the agarose solu- |
| | | • | Use bot | se a oth | transfer inside er | r pipet t nds of th | o deposit a ne bed. | a small a | mount of cooled agarose t |
| | | • Wait approximately 1 minute for the agarose | | | | | | | arose to solidify. |
| | 9. | Pou a le | ur th evel : | he c I sur | cooled ag | garose s | olution in | to the be | ed. Make sure the bed is o |
| | 10. | Allo | ow t | the | e gel to c | omplete | ely solidify. | lt will k | become firm and cool to th |

agarose is



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

14

EDVO-Kit # 371

EDVO-Kit

Module II: Separation of PCR Reactions by Electrophoresis



PREPARING THE GEL FOR ELECTROPHORESIS

11. After the gel is completely solidified, carefully and slowly remove the rubber dams or tape from the gel bed.

Be especially careful not to damage or tear the gel wells when removing the rubber dams. A thin plastic knife, spatula or pipet tip can be inserted between the gel and the dams to break possible surface tension.

12. Remove the comb by slowly pulling straight up. Do this carefully and evenly to prevent tearing the sample wells.



- 13. Place the gel (on its bed) into the electrophoresis chamber, properly oriented, centered and level on the platform.
- 14. Fill the electrophoresis apparatus chamber with the appropriate amount of diluted (1x) electrophoresis buffer.
- 15. Make sure that the gel is completely submerged under buffer before proceeding to loading the samples and conducting electrophoresis.

IMPORTANT: Check with your instructor to determine if the buffer has previously been diluted. Pour the appropriate amount of 1x buffer into the electrophoresis chamber according to Table B below.

| | Table B | Electrophoresis (Chamber) Buffer | | | | | | | | | | |
|---|-------------------|----------------------------------|-------------------------------|---|-----|--|--|--|--|--|--|--|
| Ī | ED' Ma | VOTEK odel # | Total Volume Required (ml) | Dilution 50x Conc. + Distilled Buffer (ml) + Water (m | | | | | | | | |
| | I | M6+ | 300 | 6 | 294 | | | | | | | |
| | l | MI2 | 400 | 8 | 392 | | | | | | | |
| | | M36 | 1000 | 20 | 980 | | | | | | | |



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved

For DNA analysis, the recommended electrophoresis buffer is Tris-acetate-EDTA, pH 7.8. The formula for diluting EDVOTEK (50x) concentrated buffer is one volume of buffer concentrate to every 49 volumes of distilled or deionized water. Prepare buffer as required for your electrophoresis unit. The

Experiment

Module II: Separation of PCR Reactions by Electrophoresis

Reminder:

EDVO-Kit #

37

During electrophoresis, the DNA samples migrate through the agarose gel towards the positive electrode. Before loading the samples, make sure the gel is properly oriented in the apparatus chamber.



Conducting Electrophoresis

LOAD THE SAMPLES

This experiment requires a 1.0% agarose gel. More than one group can share each gel.

Heat the 200 bp DNA ladder and 1. PCR samples for two minutes at 50°C. Allow the samples to cool for a few minutes.



- Load the DNA ladder in lane 1 of each gel. 2.
- 3. Load the entire volume (30 µl) of each PCR sample in consecutive wells.

Remember to note the wells in which your group's samples are loaded.

RUNNING THE GEL

4. After the DNA samples are loaded, carefully snap the cover down onto the electrode terminals.

Make sure that the negative and positive color-coded indicators on the cover and apparatus chamber are properly oriented.

5. Insert the plug of the black wire into the black input of the power source (negative input). Insert the plug of the red wire into the red input of the power source (positive input).

| Tab | ole C Time ar (1.0% - 7 | nd Voltage x 14 cm gel) |
|-------|----------------------------|------------------------------|
| Volts | Recomme Minimum | n ded Time Maximum |
| 125 | 55 min | l hr 15 min |
| 70 | 2 hrs 15 min | 3 hrs |
| 50 | 3 hrs 25 min | 5 hrs |

- 6. Set the power source at the required voltage and conduct electrophoresis for the length of time determined by your instructor. General guidelines are presented in Table C.
- 7. Check to see that current is flowing properly - you should see bubbles forming on the two platinum electrodes.
- 8. After the electrophoresis is completed, turn off the power, unplug the power source, disconnect the leads and remove the cover.
- Remove the gel from the bed for staining. 9.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

The Experiment

EDVO-Kit

Staining and Visualization of DNA

After electrophoresis, the agarose gels require staining in order to visualize the separated DNA samples. This experiment features a proprietary stain called InstaStain®.

INSTASTAIN® ETHIDIUM BROMIDE

Optimal visualization of PCR products on gels of 1.0% or higher concentration is obtained by staining with InstaStain® Ethidium Bromide (InstaStain® EtBr) cards.

Caution: Ethidium Bromide is a listed mutagen. Disposal of the InstaStain® EtBr cards, which contain only a few micrograms of ethidium bromide, is minimal compared to the large volume of liquid waste generated by traditional ethidium bromide staining procedures. Disposal of InstaStain® cards and gels should follow institutional guidelines for chemical waste.

INSTASTAIN® BLUE

Alternatively, InstaStain® Blue cards can be used for staining gels in this experiment. However, InstaStain® Blue is less sensitive than InstaStain® EtBr and will yield variable results.

Two options are provided for using the InstaStain® Blue cards.

Method 1: One-step Staining and Destaining with InstaStain® Blue

Method 2: Staining with InstaStain® Blue

Using Method 1, agarose gels can be stained and destained in one easy step, which can be completed in approximately 3 hours, or can be left in liquid overnight. Method 2, using InstaStain® Blue cards, requires approximately 5-10 minutes for staining. DNA bands will become visible after destaining for approximately 20 minutes, and will become sharper with additional destaining. For the best photographic results, allow the gel to destain for several hours to overnight. This will allow the stained gel to "equilibrate" in the destaining solution, resulting in dark blue DNA bands contrasting against a uniformly light blue background.

InstaStain is a registered trademark of EDVOTEK, Inc. Patents Pending.

Gels stained with InstaStain® Blue may be stored in the refrigerator for several weeks. Place the gel in a sealable plastic bag with destaining liquid. DO NOT FREEZE AGAROSE GELS! Used InstaStain® Blue cards and destained gels can be discarded in solid waste disposal. Destaining solutions can be disposed down the drain.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM The Experiment

Module II: Staining of DNA with InstaStain® Ethidium Bromide

Wear gloves and safety goggles Do not stain gel(s) in

the electrophoresis apparatus.



Visit our web site for an animated demonstration of InstaStain® EtBr. After electrophoresis, place the gel on a piece of plastic wrap on a flat surface. Moisten the gel with a few drops of electrophoresis buffer.

1

Moisten the gel.

Place the InstaStain®

Place a small weight to ensure good contact.

View on U.V. (300 nm) transilluminator

Press firmly.

card on the gel.

3

- Wearing gloves, remove the clear plastic protective sheet, and place the unprinted side of the InstaStain® EtBr card on the gel.
- Firmly run your fingers over the entire surface of the InstaStain® EtBr. Do this several times.
- 4. Place the gel casting tray and a small empty beaker on top to ensure that the InstaStain® card maintains direct contact with the gel surface.

Allow the InstaStain® EtBr card to stain the gel for 10-15 minutes.

5. After 10-15 minutes, remove the InstaStain® EtBr card. Transfer the gel to a ultraviolet (300 nm) transilluminator for viewing. Be sure to wear UV protective goggles.

DISPOSAL OF INSTASTAIN

Disposal of InstaStain® cards and gels should follow institutional guidelines for chemical waste.

Additional Notes About Staining

- If bands appear faint, or if you are not using EDVOTEK UltraSpec-Agarose[™], gels may take longer to stain with InstaStain® EtBr. Repeat staining and increase the staining time an additional 10-15 minutes.
- Gels stained alternatively with InstaStain Blue or liquid methylene blue may fade with time. Re-stain the gel to visualize the DNA bands.
- DNA 200 bp markers should be visible after staining even if the amplified DNA samples are faint or absent. If markers are not visible, troubleshoot for problems with the electrophoretic separation.



Caution: Ethidium Bromide is a listed mutagen.

Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved

EDVO-Kit #





Staining and Visualization of DNA - InstaStain® Ethidium Bromide

PHOTODOCUMENTATION OF DNA

There are many different photodocumentation systems available, including digital systems that are interfaced directly with computers. Specific instructions will vary depending upon the type of photodocumentation system you are using.

The following guidelines are for photographing gels stained with InstaStain® Ethidium Bromide, utilizing the EDVOTEK UV photodocumentation system (Cat. # 555). It is a relatively simple photodocumentation system comprised of a UV transilluminator, a 6 inch safety camera hood, and Polaroid camera fitted with a deep yellow Tiffin 40.5 mm filter. The camera uses Polaroid 667 Black and White film. The recommended settings can be used as a starting point, although optimal conditions for your system may vary.

PHOTOGRAPHY GUIDELINES

- 1. To assemble the camera, screw the handle into the center hole at the base of the camera.
- 2. Align the hood onto the camera lens.
- 3. Carefully and firmly push down the buttons on the inside of the hood on both sides of the lens.
- 4. Load your Polaroid camera with Polaroid 667 Black and White film.
- 5. Open the safety cover of the transilluminator and place the gel on the surface of the filter.
- 6. Cover the gel with the camera hood so that the hood is aligned with the camera mounting plate.
- 7. Turn on the transilluminator and photograph.
 - Recommended camera setting is f 5.6 for 2 seconds.
 - If the photograph is too light, change the aperture to f 8 and expose for 2 seconds.
 - If too dark, reduce the shutter speed to 1 second at f 5.6.

For additional information, refer to the instructions which accompany your photodocumentation system.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM The

Experiment

Staining and Visualization of DNA - InstaStain® Blue



Wear gloves and safety goggles

Do not stain gel(s) in the electrophoresis apparatus.

METHOD 1: ONE-STEP STAINING AND DESTAINING WITH INSTASTAIN® BLUE

Agarose gels can be stained and destained in one easy step with InstaStain® Blue cards. This one-step method can be completed in approximately 3 hours, or can be left overnight.

 Remove the 7 x 7 cm agarose gel from its bed and completely submerse the gel in a small, clean tray containing 75 ml of distilled or deionized water, or used electrophoresis buffer. The agarose gel should be completely covered with liquid.

Examples of small trays include large weigh boats, or small plastic food containers



- 2. Gently float a 7 x 7 cm card of InstaStain® Blue with the stain side (blue) facing the liquid.
- 3. Let the gel soak undisturbed in the liquid for approximately 3 hours. The gel can be left in the liquid overnight (cover with plastic wrap to prevent evaporation).
- 4. After staining and destaining, the gel is ready for visualization and photography.

STORAGE AND DISPOSAL OF INSTASTAIN® BLUE CARDS AND GELS

• Stained gels may be stored in the refrigerator for several weeks. Place the gel in a sealable plastic bag with destaining liquid.

DO NOT FREEZE AGAROSE GELS!

- Used InstaStain® cards and destained gels can be discarded in solid waste disposal.
- Destaining solutions can be disposed down the drain.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

20

EDVO-Kit #

Staining and Visualization of DNA - InstaStain® Blue



Wear gloves and safety goggles

METHOD 2: STAINING WITH **INSTASTAIN® BLUE CARDS**

- 1. After electrophoresis, place the agarose gel on a flat surface covered with plastic wrap.
- Wearing gloves, place the blue dye side of 2. the InstaStain® Blue card on the gel.
- Firmly run your fingers several times over 3. the entire surface of the InstaStain® card to establish good contact between the InstaStain® card and the gel.
- 4. To ensure continuous contact between the gel and the InstaStain® card, place a gel casting tray and weight, such as a small empty beaker, on top of the InstaStain® card.
- Allow the InstaStain® Blue to sit on the gel 5. for 5 to 10 minutes.
- 6. After staining, remove the InstaStain® card.

If the color of the gel appears very light, wet the gel surface with buffer or distilled water and place the InstaStain® card back on the gel for an additional 5 minutes.

Destaining and Visualization of DNA

- Transfer the gel to a large weigh boat or 7. small plastic container.
- Destain with distilled water.* 8.
 - Add approximately 100 ml of distilled water to cover the gel.

InstaStain is a registered trademark of EDVOTEK, Inc. Patents Pending.



tory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

The Biotechnology Education Company® • 1-800-EDVOTEK • www.edvotek.com

EDVO-Kit

Staining and Visualization of DNA - InstaStain® Blue

9. Repeat destaining by changing the distilled water as needed.

The larger DNA bands will initially be visible as dark blue bands against a lighter blue background. When the gel is completely destained, the larger DNA bands will become sharper and the smaller bands will be visible. With additional destaining, the entire background will become uniformly light blue.

- 10. Carefully remove the gel from the destain solution and examine the gel on a Visible Light Gel Visualization System. To optimize visibility, use the amber filter provided with EDVOTEK equipment.
- 11. If the gel is too light and bands are difficult to see, repeat the staining and destaining procedures.

* Destaining Notes

- Warmed distilled water at 37°C will accelerate destaining. Destaining will take longer with room temperature water.
- DO NOT EXCEED 37°C ! Warmer temperatures will soften the gel and may cause it to break.
- The volume of distilled water for destaining depends upon the size of the tray. Use the smallest tray available that will accommodate the gel. The gel should be completely submerged during destaining.
- Do not exceed 3 changes of water for destaining. Excessive destaining will cause the bands to be very light.

STORAGE AND DISPOSAL OF INSTASTAIN® BLUE CARDS AND GELS

• Stained gels may be stored in the refrigerator for several weeks. Place the gel in a sealable plastic bag with destaining liquid.

DO NOT FREEZE AGAROSE GELS!

- Used InstaStain® cards and destained gels can be discarded in solid waste disposal.
- Destaining solutions can be disposed down the drain.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

The Experiment



Study Questions

Answer the following study questions in your laboratory notebook or on a separate worksheet.

- 1. What is polymorphic DNA? How is it used for identification purposes?
- 2. What is CODIS? How is it used to solve crimes?
- 3. What is an STR? A VNTR? Which (STR or VNTR) is predominantly used in law enforcement? Why?





Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM



Experiment Notes



EDVOTEK - The Biotechnology Education Company® • 1-800-EDVOTEK • www.edvotek.com

Online Ordering

now available

Class size, length of laboratory sessions, and availability of equipment are factors which must be considered in the planning and the implementation of this experiment with your students. These guidelines can be adapted to fit your specific set of circumstances. If you do not find the answers to your questions in this section, a variety of resources are continuously being added to the EDVOTEK web site. In addition, Technical Service is available from 9:00 am to 6:00 pm, Eastern time zone. Call for help from our knowledgeable technical staff at 1-800-EDVOTEK (1-800-338-6835).

NATIONAL CONTENT AND SKILL STANDARDS

By performing this experiment, students will learn to extract chromosomal DNA, load samples and run agarose gel electrophoresis. Analysis of the experiments will provide students the means to transform an abstract concept into a concrete explanation. Please visit our website for specific content and skill standards for various experiments.

EDUCATIONAL RESOURCES

Electrophoresis Hints, Help and Frequently Asked Questions

EDVOTEK Electrophoresis Experiments are easy to perform and are designed for maximum success in the classroom setting. However, even the most experienced students and teachers occasionally encounter experimental

problems or difficulties. The EDVOTEK web site provides several suggestions and reminders for conducting electrophoresis, as well as answers to frequently asked electrophoresis questions.

Laboratory Extensions and Supplemental Activities

Laboratory extensions are easy to perform using EDVOTEK experiment kits. For laboratory extension suggestions, please check the EDVOTEK website, which is updated on a continuous basis with educational activities and resources.



1-800-EDVOTEK • www.edvotek.com FAX: (301) 340-0582 • email: edvotek@aol.com



EDVO-Kit # 371

Notes to the Instructor:

PCR EXPERIMENTAL SUCCESS GUIDELINES

Please refer to the Appendices section for a summary of important hints and reminders which will help maximize successful implementation of this experiment.

This experiment has two modules:

- i. PCR Amplification of Crime Scene & Suspect DNA
- II. Separation of PCR Reactions by Electrophoresis

APPROXIMATE TIME REQUIREMENTS

DNA Amplification

DNA Amplification (35 PCR cycles) will take about 70-90 minutes or can be processed overnight and held at 4°C.

Agarose Gel Preparation

There are several options for preparing agarose gels for the electrophoresis experiment. Your schedule will determine when to prepare the gel(s) for the experiment. Whether you choose to prepare the gel(s) or have the students do it, allow approximately 30-40 minutes for this procedure. Generally, 20 minutes of this time is required for gel solidification.

- Individual Gel Casting: Each student lab group can be responsible for casting their own individual gel prior to conducting the experiment.
- Batch Gel Preparation: A batch of agarose gel can be prepared for sharing by the class. To save time, a larger quantity of UltraSpec-Agarose can be prepared for sharing by the class. See instructions for "Batch Gel Preparation".
- Preparing Gels in Advance: Gels may be prepared ahead and stored for later use. Solidified gels can be stored <u>under</u> buffer in the refrigerator for up to 2 weeks.

Do not store gels at -20°C. Freezing will destroy the gels.

Gels that have been removed from their trays for storage, should be "anchored" back to the tray with a few drops of hot, molten agarose before placing the gels into the apparatus for electrophoresis. This will prevent the gels from sliding around in the trays and the chambers.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

Instructor's Guide

Notes to the Instructor:

Agarose Gel Electrophoresis

The approximate time for electrophoresis will vary from 55 minutes to 5 hours, depending on the power supply voltage.

Generally, the higher the voltage applied the faster the samples migrate. However, the maximum amount of voltage significantly depends upon the design of the electrophoresis apparatus and should not exceed manufacturers recommendations. Time and Voltage recommendations for EDVOTEK equipment are outlined in Table C.

| Tab | Die C Time ar (1.0% - 7 | nd Voltage x 14 cm gel) | | | |
|-------|----------------------------|----------------------------|--|--|--|
| Volts | Recomme Minimum | nded Time Maximum | | | |
| 125 | 55 min | l hr 15 min | | | |
| 70 | 2 hrs 15 min | 3 hrs | | | |
| 50 | 3 hrs 25 min | 5 hrs | | | |

OPTIONAL STOPPING POINTS

The experiment can be temporarily stopped after the completion of DNA Amplification (Module I) and later resumed. Experimental results will not be compromised if instructions are followed as noted under the heading "Optional Stopping Point" at the end of the procedural instructions.

STAINING AND VISUALIZATION OF PCR PRODUCTS AFTER AGAROSE GEL ELECTROPHORESIS

For this experiment, optimal visualization will be obtained by staining gels with InstaStain® Ethidium Bromide cards, which are included in this experiment. Staining of higher concentration agarose gels (1.0% or higher) require more care to obtain visible and clear results.

Disposal of the InstaStain® EtBr cards, which contain only a few micrograms of ethidium bromide, is minimal compared to the large volume of liquid waste generated by traditional ethidium bromide staining procedures. Disposal of InstaStain® cards and gels should follow institutional guideline for chemical waste.

Alternatively, InstaStain® Blue can be used for staining gels in this experiment. However, InstaStain® Blue is less sensitive than InstaStain® EtBr and will yield variable results.

LABORATORY NOTEBOOKS

It is highly recommended that students maintain a laboratory notebook to formulate hypotheses and to record experimental procedures and results.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM EDVO-Kit #

Pre-Lab Preparations

NOTE:

EDVO-Kit #

There is enough template DNA (Components D-G) for each group to perform a unique forensic crime scene scenario. There is enough material to perform 25 PCR reactions and 5 gels. Students can be divided into groups of five students per group and samples from each group can be run on a gel.

There are four different DNA samples provided. This can be an open-ended experiment that students can design – you can designate one sample as the crime scene for one group and a different sample as the crime scene for another group. Make sure to repeat the crime scene DNA as one of the suspect DNAs for the respective groups. Alternatively, students can design their own crime scene scenario and designate the DNAs for the crime scene and suspects.

- 1. For the open-ended experiments, one of the four samples has to be assembled in duplicate one of the reactions will be designated as the evidence collected at the scene of the crime and the second will be used as a possible suspect.
- 2. Thaw the frozen materials and immediately place on ice.
- 3. For the assigned samples for crime scene and suspects, aliquot and gather the materials outlined below (keep DNAs and primers on ice):

7 μl Crime scene DNA (choose DNA Template #1, #2, #3, or #4)
7 μl Suspect #1 DNA (DNA Template #1)
7 μl Suspect #2 DNA (DNA Template #2)
7 μl Suspect #3 DNA (DNA Template #3)
7 μl Suspect #4 DNA (DNA Template #4)
120 μl Primer mix
30 μl 200 base-pair ladder
5 PCR beads (in tubes)
50 μl 10X Gel Load Solution

Notes and Reminders:

- Accurate temperatures and cycle times are critical. A pre-run for one cycle (approx. 3 to 5 min) is recommended to check that the thermal cycler is properly programmed.
- For thermal cyclers which do not have a top heating plate, it is necessary to place a layer of wax above the PCR reactions in the microcentrifuge tubes to prevent evaporation. See Appendix entitled "Preparation and Handling PCR Samples with Wax ".
- Three water baths can be used for PCR if a thermal cycler is unavailable. The experiment will require great care and patience. Samples will require wax layers. See appendices entitled "Polymerase Chain Reaction Using Three Waterbaths" and "Handling samples with wax overlays".



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

EDVO-Kit #

Quantity Preparations for Agarose Gel Electrophoresis

To save time, the electrophoresis buffer and agarose gel solution can be prepared in larger quantities for sharing by the class. Unused diluted buffer can be used at a later time and solidified agarose gel solution can be remelted.

| Table D | B Ele | ulk Prepa | aration of resis Buffer |
|---------------|--------------------------------|----------------------------|----------------------------|
| Conc Buffe | entrated er (50x) + (ml) | Distilled Water (ml) | Total = Volume (ml) |
| | 60 | 2,940 | 3000 (3 L) |



Note: The UltraSpec-Agarose TM kit component is often labeled with the amount it contains. In many cases, the entire contents of the bottle is 3.0 grams. Please read the label carefully. If the amount of agarose is not specified or if the bottle's plastic seal has been broken, weigh the agarose to ensure you are using the correct amount.

BULK ELECTROPHORESIS BUFFER

Quantity (bulk) preparation for 3 liters of 1x electrophoresis buffer is outlined in Table D.

BATCH AGAROSE GELS (1.0%)

For quantity (batch) preparation of 1.0% agarose gels, see Table E.

- 1. Use a 500 ml (or larger) flask to prepare the diluted gel buffer
- Pour 4.0 grams of UltraSpec-Agarose™ into 400 ml of prepared buffer. Swirl to disperse clumps.
- 3. With a marking pen, indicate the level of solution volume on the outside of the flask.
- 4. Heat the agarose solution as outlined previously for individual gel preparation. The heating time will require adjustment due to the larger total volume of gel buffer solution.
- Cool the agarose solution to 60°C with swirling to promote even dissipation of heat. If evaporation has occurred, add distilled water to bring the solution up to the original volume as marked on the flask in step 3.
- 6. Dispense the required volume of cooled agarose solution for casting each gel. The volume required is dependent upon the size of the gel bed.
- 7. Allow the gel to completely solidify. It will become firm and cool to the touch after approximately 20 minutes. Then proceed with preparing the gel for electrophoresis.



60°C

Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved

EDVO-Kit # 371

Experiment Results and Analysis



Photos of Gel Results

Lane

- 1 200 bp DNA ladder
- 2 Crime Scene DNA
- 3 Suspect #1 DNA
- 4 Suspect #2 DNA
- Suspect #3 DNA
 Suspect #4 DNA

Note:

Depending on the PCR conditions used, a diffuse, small-molecular weight band, known as a "primer dimer", may be present below the 200 bp marker. This is a PCR artifact and can be ignored. Other minor bands may also appear due to nonspecific primer binding and the subsequent amplification of these sequences.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

Study Questions and Answers

1. What is polymorphic DNA? How is it used for identification purposes?

Polymorphic DNA refers to chromosomal regions that vary widely from person to person. This variation is usually in the length of a specific DNA region. By analyzing a number of these regions, one may obtain a "DNA fingerprint" of a person that is extremely unlikely to match the DNA fingerprint of any other individual. DNA fingerprinting is used for the identification of missing persons, human remains, and matching criminal suspects to crime scenes.

2. What is CODIS? How is it used to solve crimes?

CODIS is an acronym for the **Co**mbined **D**NA Index **S**ystem, a computerbased database containing DNA fingerprints. In the convicted offender database, DNA profiles of convicted felons are maintained. In the forensic database, DNA fingerprints from crime scenes are maintained.

3. What is an STR? A VNTR? Which (STR or VNTR) is predominantly used in law enforcement? Why?

An STR is an acronym for a short tandem repeat, a DNA sequence of 2-4 base pairs that is repeated variably from person to person. VNTRs, or variable number of tandem repeats, have longer repeat units of 15-70 base pairs. STRs are now preferred to VNTRs as the length of their amplified products requires less template DNA, often allowing even degraded samples to be amplified.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM EDVO-Kit #

EDVO-Kit #

Appendix: PCR Experimental Success Guidelines

EDVOTEK experiments which involve the extraction and amplification of DNA for fingerprinting are extremely relevant, exciting and stimulating classroom laboratory activities. These experiments have been performed successfully in many classrooms across the country, but do require careful execution because of the small volumes used. The following guidelines offer some important suggestions, reminders and hints for maximizing success.

DNA Extraction and Sample Preparation

CELL PREPARATION:

- 1. **Sufficient Cells:** It is critical that there are sufficient cells to obtain enough DNA that will yield positive DNA fingerprinting results. Cell sources include human, plant, *drosophila* and bacterial cells. Without enough cells, there will not be enough DNA template for the PCR reaction.
- 2. Human (Self) DNA Fingerprinting: Cells obtained from human sources, such as cheek cells, need to be harvested cautiously. Aerosol can result and cross-contamination among students can be a health hazard. Hair follicles do not pose the same problem and yield sufficient DNA required for the PCR reaction.
- 3. Hair Cells: At least four (4) hair follicles are needed. The **preferred** source is hair from eyebrows. Use only hairs containing a sheath, a barrel-shaped structure (often white in color) encircling the shaft near the base of the hair (see figure at left). Centrifuge the hair follicles to the bottom of the micro-



centrifuge tube to ensure direct contact with the reagents used in subsequent steps.

- 4. **Cheek Cells:** A white pellet must be visible after centrifuging the cell suspension obtained from cheek cell swabbing. If necessary, repeat the centrifugation to obtain a visible pellet. After removal of the supernatant, suspend the pellet in the chelating agent by repeated vortexing and pipetting up and down.
- 5. **Chelating Agent:** Chelating agent removes Mg (required by DNA-degrading nucleases and DNA polymerases). The small beads must be suspended in the buffer prior to delivery to the cells (i.e., mix the chelating agent just before you transfer it to the tube containing the cells.
- 6. **Boiling:** The boiling step for 10 minutes is required to obtain cell lysis. Boiling will not degrade the DNA and nucleases will NOT degrade DNA in the absence of Mg.
- 7. **Centrifugation:** Centrifuge the cell suspension carefully after cooling. If the pellet loosens, repeat this step. The supernatant should be clear, not cloudy, and the pellet should be solid at the bottom of the tube. Repeat centrifugation for a longer period of time, if necessary.
- 8. **DNA Transfer:** Transfer the DNA to a new microcentrifuge tube very carefully. It is the step prior to the PCR reaction. If any chelating agent beads (as few as one or two) are transferred, they can easily trap the Mg required by the *Taq* DNA polymerase as a cofactor for catalysis. As an additional precaution, centrifuge the supernatant a second time.

Remember: Any carry-over of chelating agent to the PCR reaction will not yield results.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM



Appendix: PCR Experimental Success Guidelines

THE PCR REACTION

- 9. Add Primers and DNA to the PCR Reaction Bead: Add the primer mixture (forward and reverse primers) and the cell DNA (supernatant) as specified in the experimental procedures to the microcentrifuge tube containing the PCR reaction bead. Make sure that the bead (which contains the *Taq* DNA polymerase, the 4XdTPs, Mg and the PCR reaction buffer) is completely dissolved. Do a quick spin in a microcentrifuge to bring the entire sample to the bottom of the tube. Prepare the control reaction similarly.
- 10. The Thermal cycler: The thermal cycler must be programmed for the correct cycle sequence. It is critical that the temperatures and the time for each of the cycles are accurate.
- 11. **Oil or Wax:** For certain thermal cyclers which do not have a top heating plate, it is necessary to overlay the reaction in the microcentrifuge tubes with oil or wax to prevent evaporation.
- 12. Manual Water Bath PCR: Three water baths can be used as an alternative to using a thermal cycler for PCR. Samples require oil or wax layers. This method requires extra care and patience and results are more variable than when using a thermal cycler.

GEL PREPARATION AND STAINING

- 13. **Concentrated agarose:** Gels of higher concentration (> 0.8%) require special attention when dissolving or re-melting. Make sure that the solution is completely clear of "clumps" or glassy granules. Distorted electrophoresis DNA band patterns will result if the gel is not properly prepared.
- 14. **Electrophoretic separation:** The tracking dye should travel at least 6 cm from the wells for adequate separation before staining.
- 15. **Staining:** Staining of higher concentration gels (> 0.8%) require additional care to obtain clear, visible results.
 - After staining (15 to 30 min.) with InstaStain® Ethidium Bromide or liquid ethidium bromide, examine the results using a UV (300nm) transilluminator. Repeat the staining as required.
 - Gels stained with InstaStain® Blue or liquid methylene blue stain may fade with time. Re-stain the gel to visualize the DNA bands.
- 16. **DNA 200 bp markers:** After staining the agarose gel, the DNA 200 bp markers should be visible after staining. If there are visible bands in the markers and control lanes, but bands in the sample lanes are faint or absent, it is possible that the DNA was not successfully extracted from the cells. If markers, control and DNA bands are all faint or absent, problems could potentially be due to improper preparation of the gel, absence of buffer in the gel, improper gel staining or a dysfunctional electrophoresis unit or power source.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved



Appendix: Polymerase Chain Reaction Using Three Waterbaths

Superior PCR results are obtained using an automated thermal cycler. However, if you do not have a thermal cycler, this experiment can be adapted to use three waterbaths (Cat. # 544). Much more care needs to be taken when using the three-waterbath PCR method. The PCR incubation sample is small and can easily be evaporated. Results using three waterbaths are often variable.

PREPARATION OF THE PCR REACTION:

- 1. The PCR reaction sample should be prepared as specified in the experiment instructions. Each PCR reaction sample contains the following three critical components:
 - PCR Reaction pellet[™]
 - (Each pellet contains *Taq* DNA polymerase, four deoxytriphosphates, Mg⁺² and buffer.)
 - Primer mix
 - DNA for amplification
- 2. After adding the components of the PCR reaction sample, use clean forceps to transfer one wax bead to the PCR tube. At the start of the PCR reaction, the wax will melt and overlay the samples to prevent evaporation during heating.

POLYMERASE CHAIN REACTION CYCLING

- 3. The three-waterbath PCR method requires three separate waterbaths, each set at different temperatures. The PCR reaction sample is sequentially cycled between the three waterbaths for a specified period of time. The sequential placement of the reaction sample in the waterbaths maintained at three different temperatures constitutes one PCR cycle. A typical PCR cycle might be set up as follows:
 - 94°C for 1 minute 45°C for 1 minute 72°C for 1 minute

It is imperative that the temperatures are accurately maintained throughout the experiment.

- 4. The PCR tube must be handled carefully when sequentially cycled between the three waterbaths. For each cycle:
 - Carefully place the PCR tube in a waterbath float. Make sure that the sample volume is at the bottom of the tube and remains undisturbed. If a tube falls on the lab bench or floor, pulse spin the tube in a balanced microcentrifuge, or shake the tube to get all of the sample to the bottom of the tube.
 - Use forceps to carefully lower the waterbath float (with tubes) sequentially into the waterbaths.
- 5. Process the PCR reaction sample for the total number of cycles specified in the experiment instructions. On the final cycle the 72°C incubation can be extended to 5 minutes.
- 6. After all the cycles are completed, the PCR sample is prepared for electrophoresis.

Please refer to the Appendix entitled "PCR Samples with Wax Overlays" for sample handling and preparation tips.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM



Appendix: Preparation and Handling of PCR Samples With Wax

For Thermal Cyclers without Heated Lids, or PCR Using Three Waterbaths

Automated thermal cyclers with heated lids are designed to surround the entire sample tube at the appropriate temperature during PCR cycles. Heating the top of the tubes during these cycles prevents the very small sample volumes from evaporating. For thermal cyclers without heated lids, or when conducting PCR by the three-waterbath method, it is necessary to add a wax bead or pellet to the reaction sample. During the PCR process, the wax will melt and overlay the samples to prevent evaporation during heating.

PREPARING THE PCR REACTION:

1. The PCR reaction sample should be prepared as specified in the experiment instructions. Each PCR reaction sample contains the following three critical components:

• PCR Reaction pellet[™] (Each pellet contains Taq DNA polymerase, four deoxytriphosphates, Mg⁺² and buffer.)

- Primer mix
- DNA for amplification
- 2. After adding the components of the PCR reaction sample, use clean forceps to transfer one wax bead to the PCR tube.
- 3. Process the PCR reaction sample for the total number of cycles specified in the experiment instructions.

PREPARING THE PCR REACTION FOR ELECTROPHORESIS:

- 4. After the cycles are completed, transfer the PCR tube to a rack and prepare the PCR sample for electrophoresis.
 - Place the PCR tube in a 94°C waterbath long enough to melt the wax overlay. Use a clean pipet to remove most of the melted wax overlay.
 - Allow a thin layer of the wax to solidify.
 - Use a clean pipet tip to gently poke a hole through the solidified wax. Remove the tip.
 - Use another clean pipet tip to enter the hole to remove the volume of mixture specified in the experiment instructions. Transfer this volume to a clean tube.
 - Add other reagents according to experiment instructions, if applicable,.
 - Add 5 µl of 10x Gel Loading solution to the sample and store on ice.
- 5. Proceed to delivery of the sample onto an agarose gel for electrophoresis as specified in the experiment instructions.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 2006, 2007, 2010, EDVOTEK, Inc., all rights reserved EVT 101121AM

| EDVOTEK. | May be used Standard. 29 | lateria I to com CFR 19 s | al Safety Data Si bly with OSHA's Hazard 10.1200 Standard mus pecific requirements. | heet I Communicatio st be consulted t | n for |
|--|-----------------------------|------------------------------------|--|---|-----------------|
| IDENTITY (As Used on Label and List) | | | Note: Blank spaces are | not permitted. If | any item is not |
| Agarose | | | be marked to indicate th | ation is available, iat. | the space must |
| Section I | | | | | |
| Manufacturer's Name | | Emer | gency Telephone Nur | nber (301) 2 | 51-5990 |
| EDVOTEK, Inc. | | Talan | ana Number for inform | otion | |
| Address (Number, Street, City, State, | Zip Code) | Telepi | | (301) 2 | 51-5990 |
| 14676 Rothgeb Drive | | Date | 11/21/10 11/21/10 | 0 | |
| Rockville, MD 20850 | | Signa | ture of Preparer (option | al) | |
| Section II - Hazardous Ingred | lients/Iden | l tify Ir | formation | | |
| Chemical Identity; Common Name(s)] | OSHA | PEL | O ACGIH TLV Re | commended | % (Optional) |
| Standard. CAS #9012-36-6 Section III - Physical/Chemic | al Charact | teristi | | | n |
| Boiling Point For 1% solution | 94 F | Spe | cific Gravity ($H_2^0 = 1$) | | No data |
| Vapor Pressure (mm Hg.) | No data | Mel | Melting Point | | No data |
| Vapor Density (AIR = 1) | No data | Eva (Bu | poration Rate tyl Acetate = 1) | No data | |
| Solubility in Water Insoluble - cold | 1 | | | | |
| Appearance and Odor White p | owder, no odo | r | | | |
| Section IV - Physical/Chemic | al Charac | terist | cs N.D. = No da | ta | |
| Flash Point (Method Used) No data | ì | Flan | mable Limits | LEL N.D. | UEL N.D. |
| Extinguishing Media Water spray, dr | y chemical, ca | rbon di | oxide, halon or standa | ard foam | |
| Special Fire Fighting Procedures Possible fire ha | zard when ex | posed to | heat or flame | | |
| Unusual Fire and Explosion Hazards | None | | | | |

| Unstable Stable available rproducts May Occur Will Not Occur azard Data | X | Conditi | ons to Avoid None | |
|---|---|---|---|--|
| Stable available rproducts May Occur Will Not Occur azard Data | X | | None | |
| available /products May Occur Will Not Occur azard Data | | | | |
| May Occur Will Not Occur Azard Data | | | | |
| May Occur Will Not Occur azard Data | | | | |
| Will Not Occur | | Condit | tions to Avoid | |
| azard Data | Х | 1 | None | |
| | | | | |
| Inhalatio | 1? Yes | | Skin? Yes | Ingestion? |
| hronic) No data availabi | e Iı | ngestion | 1: Large amounts | may cause diarrhea |
| NTP? | 11 | IAF | C Monographs? | OSHA Regulation? |
| osure No data | availabl | le | | |
| y Aggravated by | Exposu | re No | data available | |
| | | 110 | und a value of the | |
| Treat syr | nptomat | ically a | nd supportively | |
| ons for Safe | Hand | ling a | nd Use | |
| aterial is Release | ed for Sp | oilled | | |
| Sweep up and J | place in s | suitable | container for disp | osal |
| | | | | |
| Normal solid w | aste disp | posal | | |
| andling and Stor | ing | | | |
| None | | | | |
| | | | | |
| None | | | | |
| Measures | | | | |
| cify Type) Che | mical ca | rtridge | respirator with ful | l facepiece. |
| Local Exhaust | | | Specia | al |
| Mechanical (Ger | neral)Gei | n. diluti | on ventilationOthe | r |
| | | | Eye Protection | Splash proof goggles |
| Equipment Ir | nperviou | is clothi | ing to prevent skin | contact |
| N | one | | | |
| | hronic) No data availabi NTP? Isure No data (Aggravated by ures Treat syr Das for Safe aterial is Release Sweep up and p Normal solid w andling and Stor None Measures Sify Type) Che Local Exhaust Mechanical (Ger Equipment Ir N | hronic) No data available II NTP? Sure No data available / Aggravated by Exposu ares Treat symptomat Treat symptomat Das for Safe Hand aterial is Released for Sp Sweep up and place in si Normal solid waste disp andling and Storing None None Measures Sify Type) Chemical ce Local Exhaust Mechanical (General)Ge Equipment Imperviou None | hronic) No data available Ingestion NTP? IAF Ssure No data available (Aggravated by Exposure No ures Treat symptomatically a Dns for Safe Handling a atterial is Released for Spilled Sweep up and place in suitable Normal solid waste disposal andling and Storing None None Measures Sify Type) Chemical cartridge Local Exhaust Mechanical (General)Gen. diluti Equipment Impervious cloth None | Tes Tes Nornic) No data available Ingestion: Large amounts in NTP? IARC Monographs? Isure No data available / Aggravated by Exposure No data available / Aggravated by Exposure No data available Treat symptomatically and supportively Dns for Safe Handling and Use atterial is Released for Spilled Sweep up and place in suitable container for disp Normal solid waste disposal andling and Storing None Measures iffy Type) Chemical cartridge respirator with full Local Exhaust Specia Mechanical (General)Gen. dilution ventilation for the Eye Protection Equipment Impervious clothing to prevent skin |

| ED VOTEK. | May be used Standard. 29 | to com CFR 19 s | al Safety Data S bly with OSHA's Hazard 10.1200 Standard mus pecific requirements. | heet d Communicatio st be consulted t | n for |
|---|-----------------------------|-----------------------|---|--|-----------------------------------|
| IDENTITY (As Used on Label and List) Tris Buffer | | | Note: Blank spaces are applicable, or no inform be marked to indicate th | e not permitted. If ation is available, nat. | any item is not the space must |
| Section I | | | | | |
| Manufacturer's Name | | Emer | gency Telephone Nu | mber (301) 2 | 51-5990 |
| EDVOTEK, Inc. | | Tolopi | ana Number for inform | (001) 2 | .01-0000 |
| Address (Number, Street, City, State, | Zip Code) | Telepi | ione Number for Inform | (301) 2 | 51-5990 |
| 14676 Rothgeb Drive | | Date | Prepared | 1/10 | |
| Rockville, MD 20850 | | Signa | ture of Proparer (option | | |
| | | Signa | ure of r reparer (option | ai) | |
| Section II - Hazardous Ingred | lients/Iden | tifv Ir | formation | | |
| Hazardous Components [Specific | 0011 | | 0 | ther Limits | 0((Ontine - 1) |
| Tris (hydroxymethyl) aminoethane | OSHA | PEL | ACGIH ILV Re | commended | % (Optional) |
| CAS# 77-86-1 | | | | | |
| | | | | | |
| Section III - Physical/Chemic | al Charact | eristi | cs | | |
| Boiling Point | No data | Spe | cific Gravity ($H_20 = 1$) |) | Not availab |
| Vapor Pressure (mm Hg.) | Negligible | Mel | ting Point | | 171°C |
| Vapor Density (AIR = 1) | No data | Eva (Bu | Evaporation Rate (Butyl Acetate = 1) | | Negligible |
| Solubility in Water Sol | uble | | | | |
| Appearance and Odor Wh | ite crystalline | e solid | | | |
| Section IV - Physical/Chemic | al Charac | teristi | cs | | |
| Flash Point (Method Used) | | Flam | mable Limits | LEL | UEL |
| Extinguishing Media Water spray | /, dry chemic | al, car | oon dioxide | • | |
| Special Fire Fighting Procedures | | | | | |
| Wear SCBA | A apparatus a | and pro | tective clothing | | |
| Unusual Fire and Explosion Hazards | | | | | |
| Fire or exce | essive heat n | nay pro | ducehazardous de | composition p | products. |

| Section V - React | tivity Data | | | | | |
|--------------------------------|---|-------------|-------------|-----------------|---------------|------------------|
| Stability | Unstable | | Conditi | ons to Avoid | | |
| | Stable | Х | 1 | | | |
| Incompatibility | | Stron | ıg oxidize | rs | | |
| Hazardous Decompositio Co | on or Byproducts ombustion will prod | uce carbo | on dioxide | and probably | y carbon mo | noxide. |
| Hazardous | May Occur | | Condi | ions to Avoid | | |
| Polymerization | Will Not Occu | ır X | | | | |
| Section VI - Heal | th Hazard Data | a | | | | |
| Route(s) of Entry: | Inhala Yes | ation? | | Skin? Yes | | Ingestion? |
| Health Hazards (Acute | and Chronic) | ow harard | l/inhalatio | n, skin, inges | stion | |
| Carcinogenicity: | NTF | ?? | IAF | C Monograp | ohs? | OSHA Regulation? |
| Signs and Symptoms | of Exposure | | | No data | | |
| Medical Conditions Ge | enerally Aggravated | by Expos | sure | None | | |
| Emergency First Aid F | Procedures | | | | | |
| Treat | | | | | | |
| 11041 | symptomatically a | ia suppor | tivery | | | |
| Section VII - Pred | cautions for Sa | ife Han | dling a | nd Use | | |
| Steps to be Taken in c | ase Material is Rele | ased for | Spilled | | | |
| Venti | late area and wash s | spill area | | | | |
| Waste Disposal Metho Dispo | od ose by incineration of | or with lie | censes che | emical waste | disposal | |
| | | | | | | |
| Precautions to be Take None | n in Handling and S | Storing | | | | |
| Other Precautions None | | | | | | |
| Section VIII - Cor | ntrol Measures | ; | | | | |
| Respiratory Protection | n (Specify Type) | Supplim | ientary ve | ntilation or re | espiratory pr | rotection |
| Ventilation | Local Exhaus | st | Yes | | Special | |
| Ventilation | Mech | anical (Gr | eneral) | Yes | Other | |
| Protective Gloves | Yes | un | 5110.2.5, | Eye Protec | tion | Yes |
| Other Protective Cloth | ing or Equipment | None | | | | |
| Work/Hygienic Practic | es | None | | | | |

| | | | | | Section V - Reactivit | y Data | | | | |
|--|--|--|---|---|--|--|--|---|--|--|
| | May be used | aterial Safety Data S to comply with OSHA's Hazan | heet | n | Stability | Unstable | Cond | itions to Avoid | | |
| EDVOTEK. | Standard. 29 | CFR 1910.1200 Standard mu specific requirements. | st be consulted f | for | Incompatibility | Stable Strong oxidizin | x g agents | None | | |
| IDENTITY (As Used on Label and List) | | Note: Blank spaces are | e not permitted. If | any item is not | Hazardous Decomposition or | Byproducts Carb | on monoxide | , Carbon dio | xide | |
| 50x Electrophoresis | Buffer | be marked to indicate the | hat. | the space must | Hazardous | May Occur | Con | ditions to Avoid | d | |
| Section I | | Emergen av Telenhause Niv | | | Polymerization | Will Not Occur | Х | None | | |
| Manufacturer's Name | | Emergency relephone Nu | (301) 2 | 51-5990 | Section VI - Health H | Hazard Data | | | | |
| Address (Number, Street, City, State, | Zip Code) | Telephone Number for inform | nation (301) 2 | 51-5990 | Route(s) of Entry: | Inhalation' | ? Yes | Skin? Y | es | Ingestiop? |
| 14676 Bothgeb Drive | | Date Prepared | 11/21/10 | | Health Hazards (Acute and | None | | | | |
| Rockville, MD 20850 | | Signature of Preparer (option | nal) | | Carcinogenicity: None ider | ntified NTP? | 1A | ARC Monogra | aphs? | OSHA Regulation? |
| | | | | | Signs and Symptoms of Ex | posure Irritation | to upper resp | piratory tract, | , skin, eyes | |
| Section II - Hazardous Ingred | lients/Iden | tify Information | | | Medical Conditions Genera | ally Aggravated by E | xposure | None | | |
| Hazardous Components [Specific Other Limits Chemical Identity; Common Name(s)] OSHA PEL ACGIH TLV Recommended % (Optional) | | % (Optional) | Francisco First Aid David | | 70 | None | | | | |
| This product contains no hazardo Communication Standard. | ous materials a | is defined by the OSHA Haz | ard | | Eyes: Flush with water | Inhalation: Move | e to fresh air | s, give large Skin: W | amounts of wa Vash with soap | and water |
| | | | | | Section VII - Precaut | tions for Safe I | Handling | and Use | | |
| Section III - Physical/Chemic | al Charact | eristics | | | Steps to be Taken in case N | Material is Released | for Spilled | Wear suitab | le protective c | lothing. Mop up spill |
| Boiling Point | No data | Specific Gravity (H ₂ 0 = 1 |) | No data | and rinse w | ith water, or collect | in absorptive | e material and | d dispose of the | e absorptive material. |
| Vapor Pressure (mm Hg.) | No data | Melting Point | | No data | Waste Disposal Method | Dispose in accorda enviromental regul | nce with all ations. | applicable fe | deral, state, an | .d local |
| Vapor Density (AIR = 1) | No data | Evaporation Rate (Butyl Acetate = 1) | | No data | Precautions to be Taken in | Handling and Storir | ng | | | |
| Solubility in Water Appreciable, (| greater than 10 | 0%) | | | | Avoid eye and skir | n contact. | | | |
| Appearance and Odor Clear, liquid, s | light vinegar o | odor | | | Other Precautions | N. | | | | |
| Section IV - Physical/Chemic | al Charact | teristics N.D. = N | No data | | Section VIII - Contro | None | | | | |
| Flash Point (Method Used) | | Flammable Limits | LEL | UEL | Respiratory Protection (Sp | ecify Type) | | | | |
| Extinguishing Media | | | N.D. | N.D. | Ventilation | Local Exhaust | Yes | | Special N | Jone |
| | Jse extinguish | ing media appropriate for su | rrounding fire. | | | Mechanical (Gene | eral) Yes | _ | Other N | √one |
| Special Fire Fighting Procedures V | Vear protective | e equipment and SCBA with | full facepiece | | Protective Gloves Yes | | | Eye Prote | ction Safety g | goggles |
| | perated in pos | itive pressure mode. | | | Other Protective Clothing o | r Equipment No | me | | | |
| Unusual Fire and Explosion Hazards | Jone identified | 1 | | | Work/Hygienic Practices | No | | | | |
| | | | | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | INC | lie | | | |
| | | | | | Section V - Reactivity | y Data | | | | |
| | May be used | laterial Safety Data S | heet | | Stability | Unstable | Condi | tions to Avoid | | |
| EDVOTEK. | Standard. 29 | CFR 1910.1200 Standard mu | st be consulted f | n for | | Stable X | | None | | |
| | | specific requirements. | | | Incompatibility | None known | | | | |
| IDENTITY (As Used on Label and List) | ncentrate 1(| Note: Blank spaces an applicable, or no inform | e not permitted. If nation is available, | any item is not the space must | Hazardous Decomposition or E | Byproducts Sulfur oxides and | l bromides | | | |
| Section I | neentrate, re | be marked to indicate t | nat. | | Hazardous | May Occur | Cond | itions to Avoid | | |
| Manufacturer's Name | | | | | Polymorization | | Λ | None | | |
| | | Emergency Telephone Nu | mber (201) 2 | 51-5000 | Polymerization | Will Not Occur | | | | |
| EDVOTEK, Inc. | Zin Code) | Emergency Telephone Nu Telephone Number for inform | mber (301) 2 | 251-5990 | Section VI - Health H Route(s) of Entry: | Will Not Occur lazard Data Inhalation? | | Skin? | Vac | Ingestion? |
| EDVOTEK, Inc. Address (Number, Street, City, State, | Zip Code) | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 | mber (301) 2 nation (301) 2 | 251-5990 | Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eve contact: May | Will Not Occur lazard Data Inhalation? Yes Chronic) | No data avai | Skin? Yes lable for othe | Yes | Ingestion? |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 | Zip Code) | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (optior | (301) 2 (301) 2 (301) 2 (301) 2 | 251-5990 | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation NTP? No data | No data avaii IA No da | Skin? Yes lable for othe RC Monogra | Yes r routes phs? (| Ingestion? |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 | Zip Code) | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option | mber (301) 2 nation (301) 2 | 251-5990 | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exp Max cause skin or eve ir | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation NTP? No data posure ritation | No data avai IA No da | Skin? Yes lable for othe RC Monograp | Yes r routes phs? (No | Ingestion? DSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingreet Hazardous Components (Specific | Zip Code) | Emergency Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information | mber (301) 2 nation (301) 2 nal) | 251-5990 251-5990 | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exp May cause skin or eye in Medical Conditions Genera None reported | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation NTP? No data posure rritation Illy Aggravated by E | No data avail IA No da xposure | Skin? Yes lable for othe RC Monogra _i tta | Yes r routes phs? (No | Ingestion? DSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingree Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard | Zip Code) dients/Iden OSHA ous material | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information APEL ACGIH TLV Re s as defined by the OSH | mber (301) 2 nation (301) 2 nal) nal) Dther Limits acommended A Hazard | 251-5990 251-5990 % (Optional) | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exp May cause skin or eye in Medical Conditions Genera None reported Emergency First Aid Procee | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation NTP? No data posure rritation illy Aggravated by E dures | No data avaii IA No da xposure | Skin? Yes lable for othe RC Monogra _i tta | Yes r routes phs? (No | Ingestion? OSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingred Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. | Zip Code) dients/Iden OSHA ous material | Emergency Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information NPEL ACGIH TLV Res as defined by the OSH | mber (301) 2 nation (301) 2 nal) Dther Limits acommended IA Hazard | 251-5990 251-5990 % (Optional) | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Proce Treat symp Rinse cont | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation mTP? No data posure rritation IIIy Aggravated by E dures ptomatically and sup tacted area with cop | No data avail IA No da xposure poortively ious amounts | Skin? Yes lable for othe RC Monograp ata | Yes r routes phs? (No | Ingestion? DSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingree Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic | Zip Code) dients/Iden OSHA ous material | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information APEL ACGIH TLV Re Is as defined by the OSH teristics | mber (301) 2 nation (301) 2 nal) Dther Limits commended IA Hazard | 251-5990 (51-5990 % (Optional) | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Procee Treat symp Rinse cont Section VII - Precaut Cutoes to be Tokme in eacob | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation NTP? No data posure rritation Illy Aggravated by E dures ptomatically and sup tacted area with cop tions for Safe | No data avail IA No da xposure poortively ious amounts landling a for Spillod | Skin? Yes lable for othe RC Monograp tta s of water. and Use | Yes r routes phs? (No | Ingestion? DSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingred Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic Boiling Point | Zip Code) dients/Iden OSHA ous material cal Charact | Emergency Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information APEL ACGIH TLV Rev as as defined by the OSH teristics Specific Gravity (H ₂ 0 = 1 | (301) 2 nation (301) 2 nal) Dther Limits accommended [A Hazard) | 251-5990 251-5990 % (Optional) No data | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exy May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Procee Treat symp Rinse cont Steps to be Taken in case M Rinse contacted area with c | Will Not Occur Hazard Data Inhalation? Yes Chronic) cause irritation No data posure rritation Ully Aggravated by E dures ptomatically and sup tacted area with cop tions for Safe H Aaterial is Released copious amounts of | No data avaii IA No da xposure opportively ious amounts tandling a for Spilled water. | Skin? Yes lable for othe RC Monograp ta | Yes r routes phs? C No | Ingestion? DSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingred Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic Boiling Point Vapor Pressure (mm Hg.) | Zip Code) dients/Iden OSHA ous material al Charact No data | Emergency Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information APEL ACGIH TLV Reference of the action of th | mber (301) 2 nation (301) 2 nal) Dther Limits ccommended IA Hazard | 251-5990 251-5990 % (Optional) No data N/A | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Procee Treat symp Rinse cont Section VII - Precaut Steps to be Taken in case M Rinse contacted area with c Waste Disposal Method Observe all Enderg | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation NTP? No data posure ritation ally Aggravated by E dures ptomatically and sup tacted area with cop tions for Safe H daterial is Released copious amounts of the multiple of the superstances the sup | No data avaii IA No da xposure portively ious amounts for Spilled water. | Skin? Yes lable for othe RC Monograp ta s of water. | Yes r routes phs? (No | Ingestion? DSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingreet Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic Boiling Point Vapor Pressure (mm Hg.) Vapor Density (AIR = 1) | Zip Code) dients/Iden OSHA ous material al Charact No data No data | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information PEL ACGIH TLV Re is as defined by the OSH seristics Specific Gravity (H ₂ 0 = 1 Melting Point Evaporation Rate (Butyl Acetate = 1) | mber (301) 2 nation (301) 2 nal) Dther Limits ccommended IA Hazard | 251-5990 251-5990 % (Optional) % (Optional) No data N/A No data | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Procee Treat symp Rinse cont Steps to be Taken in case M Rinse contacted area with co Waste Disposal Method Observe all federal, state, a Precautions to be Taken in Laken in Laken | Will Not Occur Hazard Data Inhalation? Yes Chronic) cause irritation NTP? No data posure mitation Ully Aggravated by E dures ptomatically and sup tacted area with cop tions for Safe H Aaterial is Released copious amounts of the mitation | No data avaii IA No da xposure poprtively ious amounts landling a for Spilled water. | Skin? Yes lable for othe RC Monograp ta | Yes r routes phs? C No | Ingestion? OSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingred Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemica Boiling Point Vapor Pressure (mm Hg.) Vapor Density (AIR = 1) Solubility in Water soluble | Zip Code) dients/Iden OSHA ous material al Charact No data No data No data | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information APEL ACGIH TLV Re Is as defined by the OSH teristics Specific Gravity (H20 = 1 Melting Point Evaporation Rate (Butyl Acetate = 1) | mber (301) 2 nation (301) 2 nal) Dther Limits ecommended IA Hazard | 251-5990 251-5990 % (Optional) No data N/A No data | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Procee Treat symp Rinse cont Section VII - Precaut Steps to be Taken in case M Rinse contacted area with c Uwaste Disposal Method Observe all federal, state, a Precautions to be Taken in I Avoid eye and skin contact | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation NTP? No data posure ritation lily Aggravated by E dures ptomatically and sup tacted area with cop tions for Safe H Aaterial is Released copious amounts of und local regulations Handling and Storin | No data avaii IA No da xposure portively ious amounts landling a for Spilled water. g | Skin? Yes lable for othe RC Monograp RC Monograp ta | Yes r routes phs? (No | Ingestion? |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingreet Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic Boiling Point Vapor Pressure (mm Hg.) Vapor Density (AIR = 1) Solubility in Water Soluble Appearance and Odor Blue Liquid, no o | Zip Code) dients/Iden OSHA ous material al Charact No data No data No data | Emergency Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information PEL ACGIH TLV Re is as defined by the OSH is eristics Specific Gravity (H ₂ 0 = 1 Melting Point Evaporation Rate (Butyl Acetate = 1) | mber (301) 2 nation (301) 2 nal) Dther Limits ccommended IA Hazard | 251-5990 51-5990 % (Optional) % (Optional) No data N/A No data | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Procee Treat symp Rinse cont Section VII - Precauti Steps to be Taken in case M Rinse contacted area with c Waste Disposal Method Observe all federal, state, a Precautions to be Taken in I Avoid eye and skin contact Other Precautions None | Mill Not Occur lazard Data Inhalation? Yes Chronic) cause irritation NTP? No data posure rritation Illy Aggravated by E dures ptomatically and sup tacted area with cop tions for Safe H Adaterial is Released copious amounts of the and local regulations Handling and Storin | No data avaii IA No da xposure poprtively ious amounts fandling a for Spilled water. 9 | Skin? Yes lable for othe RC Monograp ta | Yes r routes phs? C No | Ingestion? OSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingreet Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic Boiling Point Vapor Pressure (mm Hg.) Vapor Density (AIR = 1) Solubility in Water Soluble Appearance and Odor Blue liquid, no of Section IV - Physical/Chemic | Zip Code) dients/Iden OSHA ous material al Charact No data No data No data odor cal Characc | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information APEL ACGIH TLV Re is as defined by the OSH is as defined by the OSH Everistics Specific Gravity (H ₂ 0 = 1 Melting Point Evaporation Rate (Butyl Acetate = 1) teristics | mber (301) 2 nation (301) 2 nal) Dther Limits ccommended IA Hazard | 251-5990 51-5990 % (Optional) % (Optional) No data N/A No data | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Procee Treat symp Rinse cont Section VII - Precauti Steps to be Taken in case N Rinse contacted area with co Waste Disposal Method Observe all federal, state, a Precautions to be Taken in I Avoid eye and skin contact. Other Precautions None | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation NTP? No data posure rritation Illy Aggravated by E dures ptomatically and sup tacted area with cop itions for Safe H Material is Released copious amounts of the Inhalation of the selessed multiple of the selessed tacted area with cop tions for Safe H Material is Released tacted area units of the Inhalation of the selessed tacted area units of the Inhalation of the selessed tacted area units of the Material is Released tacted area units of the Material is Released tacted area units of the Inhalation of the selessed tacted area units of the selessed tacted area | No data avaii IA No da xposure opportively ious amounts tandling a for Spilled water. g | Skin? Yes lable for othe RC Monograp ta | Yes r routes phs? C No | Ingestion? DSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingreet Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic Boiling Point Vapor Pressure (mm Hg.) Vapor Density (AIR = 1) Solubility in Water Soluble Appearance and Odor Blue liquid, no co Section IV - Physical/Chemic Flash Point (Method Used) | Zip Code) dients/Iden OSHA ous material al Charact No data No data No data No data | Emergency Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information APEL ACGIH TLV Re as as defined by the OSH eristics Specific Gravity (H ₂ 0 = 1 Melting Point Evaporation Rate (Butyl Acetate = 1) teristics Flammable Limits | mber (301) 2 nation (301) 2 nal) Dther Limits accommended IA Hazard | 251-5990 51-5990 % (Optional) % (Optional) No data N/A No data | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Proceet Treat symp Rinse cont Section VII - Precauti Steps to be Taken in case N Rinse contacted area with co Usate Disposal Method Observe all federal, state, a Precautions to be Taken in I Avoid eye and skin contact. Other Precautions None Section VIII - Control Bespiratory Protection (Per | Will Not Occur Inhalation? Yes Chronic) cause irritation NTP? No data posure ritation ully Aggravated by E dures ptomatically and sup tions for Safe H Adaterial is Released copious amounts of the und local regulations Handling and Storin Handling and Storin | No data avaii IA No da xposure opportively ious amounts tandling a for Spilled water. g | Skin? Yes lable for othe RC Monograp ta s of water. and Use | Yes r routes phs? C No | Ingestion? DSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingree Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic Boiling Point Vapor Pressure (mm Hg.) Vapor Density (AIR = 1) Solubility in Water Soluble Appearance and Odor Blue liquid, no c Section IV - Physical/Chemic Flash Point (Method Used) No data Extinguishing Media | Zip Code) dients/Iden OSHA OUS material al Charact No data No data No data dor cal Charact | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information APEL ACGIH TLV Re Is as defined by the OSH teristics Specific Gravity (H20 = 1 Melting Point Evaporation Rate (Butyl Acetate = 1) teristics Flammable Limits | mber (301) 2 nation (301) 2 nal) Dther Limits commended IA Hazard)) | 251-5990 251-5990 % (Optional) % (Optional) No data N/A No data | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Procee Treat symp Rinse cont Section VII - Precaut Steps to be Taken in case M Rinse contacted area with c Usaste Disposal Method Observe all federal, state, a Precautions to be Taken in I Avoid eye and skin contact. Other Precautions None Section VIII - Control Respiratory Protection (Spe Ventilation | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation models posure ritation ally Aggravated by E dures ptomatically and sup tacted area with cop tions for Safe H Aaterial is Released copious amounts of the multiple of the second multiple | No data avaii IA No da xposure portively ious amounts tandling a for Spilled water. g | Skin? Yes lable for othe RC Monograp RC Monograp at a s of water. and Use respirator wi | Yes r routes phs? C No No Special | Ingestion? OSHA Regulation? data |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingree Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic Boiling Point Vapor Pressure (mm Hg.) Vapor Density (AIR = 1) Solubility in Water Soluble Appearance and Odor Blue liquid, no c Section IV - Physical/Chemic Flash Point (Method Used) No data Extinguishing Media Dry chemical, c | Zip Code) dients/Iden OSHA OUS material al Charact No data No data No data dor cal Characc arbon dioxid | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information APEL ACGIH TLV Re Is as defined by the OSH teristics Specific Gravity (H ₂ 0 = 1 Melting Point Evaporation Rate (Butyl Acetate = 1) teristics Flammable Limits le, water spray or foam | mber (301) 2 nation (301) 2 nal) Dther Limits commended IA Hazard)) | 251-5990 251-5990 % (Optional) % (Optional) No data N/A No data | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Procee Treat symp Rinse contacted area with c Steps to be Taken in case M Rinse contacted area with c Uses Disposal Method Observe all federal, state, a Precautions to be Taken in I Avoid eye and skin contact Other Precautions None Section VIII - Control Respiratory Protection (Spe Ventilation | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation modelships posure ritation MTP? No data posure printation dures pomatically and sup tacted area with cop tions for Safe H daterial is Released copious amounts of the modelships and local regulations Handling and Storin Handling and Storin Handling and Storin Local Exhaust Mechanical (Gene | No data avaii IA No da xposure portively ious amounts landling a for Spilled water. g | Skin? Yes lable for othe RC Monograp s of water. and Use respirator wi (cs fes | Yes r routes phs? C No No Special Other | Ingestion? OSHA Regulation? data OSHA Regulation? data OSHA Regulation? data OSHA Regulation? OSHA Regulatio |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingreet Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic Boiling Point Vapor Pressure (mm Hg.) Vapor Density (AIR = 1) Solubility in Water Solubile Appearance and Odor Blue liquid, no c Section IV - Physical/Chemic Flash Point (Method Used) No data Extinguishing Media Dry chemical, c Special Fire Fighting Procedures Use agents suita | Zip Code) dients/Iden OSHA OUS material al Charact No data No data No data odor cal Characc arbon dioxid ble for type | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information APEL ACGIH TLV Reference is as defined by the OSH teristics Specific Gravity (H20 = 1 Melting Point Evaporation Rate (Butyl Acetate = 1) teristics Flammable Limits le, water spray or foam of surrounding fire. Ket | mber (301) 2 nation (301) 2 nal) Dther Limits commended IA Hazard)) LEL No data | 251-5990 251-5990 % (Optional) % (Optional) No data N/A No data | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Proceut Rinse contacted area with c Steps to be Taken in case M Rinse contacted area with c Waste Disposal Method Observe all federal, state, a Precautions to be Taken in I Avoid eye and skin contact Other Precautions None Section VIII - Control Respiratory Protection (Spe Ventilation Protective Gloves yes | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation yosure ritation ully Aggravated by E dures pomatically and suptacted area with cop taterial is Released copious amounts of the and local regulations Handling and Storin Linadium Chronical Local Exhaust Mechanical (Geness | No data avaii IA No data xposure portively ious amounts landling a for Spilled water. 9 9 | Skin? Yes lable for othe RC Monograp s of water. and Use respirator wi (cs (cs Eye Protec | Yes r routes phs? C No No Special Other tionSplash prov | Ingestion? OSHA Regulation? data OSHA Regulation? data or cartridge. Yes None of goggles |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingreet Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic Boiling Point Vapor Pressure (mm Hg.) Vapor Density (AIR = 1) Solubility in Water Solubility in Water Solubile Appearance and Odor Blue liquid, no co Section IV - Physical/Chemic Flash Point (Method Used) No data Extinguishing Media Dry chemical, co Special Fire Fighting Procedures Use agents suita breathing hazard | Zip Code) dients/Iden OSHA OUS material al Charact No data No data No data ador cal Charact arbon dioxid ble for type dous sulfur o | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information APEL ACGIH TLV Reference is as defined by the OSH teristics Specific Gravity (H20 = 1 Melting Point Evaporation Rate (Butyl Acetate = 1) teristics Flammable Limits le, water spray or foam of surrounding fire. Ket ixides and bromides. Wo | mber (301) 2 nation (301) 2 nal) Dther Limits commended IA Hazard)) LEL No data ep upwind, av | 251-5990 251-5990 % (Optional) % (Optional) No data N/A No data UEL No data | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Proceut Rinse contacted area with c Steps to be Taken in case M Rinse contacted area with c Waste Disposal Method Observe all federal, state, a Precautions to be Taken in I Avoid eye and skin contact Other Precautions None Section VIII - Control Respiratory Protection (Spe Ventilation Protective Gloves yes Other Protective Clothing on | Will Not Occur lazard Data Inhalation? Yes Chronic) cause irritation yosure ritation ully Aggravated by E dures pomatically and suptacted area with cop taterial is Released copious amounts of the and local regulations Handling and Storin Local Exhaust Mechanical (Gene s r Equipment | No data avaii IA No data xposure portively ious amounts landling a for Spilled water. 9 9 | Skin? Yes lable for othe RC Monograp s of water. and Use respirator wi (cs (cs Eye Protec | Yes r routes phs? C No No Special Other tionSplash provi | Ingestion? DSHA Regulation? data DSHA Regulation? data DSHA Regulation? data DSHA Regulation? data DSHA Regulation? DSHA Regu |
| EDVOTEK, Inc. Address (Number, Street, City, State, 14676 Rothgeb Drive Rockville, MD 20850 Section II - Hazardous Ingreet Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazard Communication Standard. Section III - Physical/Chemic Boiling Point Vapor Pressure (mm Hg.) Vapor Density (AIR = 1) Solubility in Water Soluble Appearance and Odor Blue liquid, no co Section IV - Physical/Chemic Flash Point (Method Used) No data Extinguishing Media Dry chemical, cc Special Fire Fighting Procedures Use agents suita breathing hazard Unusual Fire and Explosion Hazards | Zip Code) dients/Iden OSHA ous material al Charact No data No data No data No data door cal Charact cal Charact ble for type lous sulfur of | Emergency Telephone Nu Telephone Number for inform Date Prepared 11/21/10 Signature of Preparer (option tify Information PEL ACGIH TLV Re- is as defined by the OSH is as defined by the OSH is eristics Specific Gravity (H ₂ 0 = 1 Melting Point Evaporation Rate (Butyl Acetate = 1) teristics Flammable Limits le, water spray or foam of surrounding fire. Kee oxides and bromides. We | mber (301) 2 nation (301) 2 nal) Dther Limits commended IA Hazard) LEL No data ep upwind, av ear SCBA. | 251-5990 251-5990 % (Optional) % (Optional) No data N/A No data UEL No data | Polymerization Section VI - Health H Route(s) of Entry: Health Hazards (Acute and Acute eye contact: May Carcinogenicity: None Signs and Symptoms of Exq May cause skin or eye ir Medical Conditions Genera None reported Emergency First Aid Proceet Treat symp Rinse cont Section VII - Precaution Vaste Disposal Method Observe all federal, state, a Precautions to be Taken in I Avoid eye and skin contact Other Precautions None Section VIII - Control Respiratory Protection (Spr Ventilation Protective Gloves yes Other Protective Clothing of Work/Hyvienic Prartinge | Will Not Occur Inhalation? Yes Chronic) cause irritation yosure rintalion? posure illy Aggravated by E dures ptomatically and suptacted area with cop ions for Safe H Adaterial is Released copious amounts of the ind local regulations Handling and Storin Local Exhaust Mechanical (Geness r Equipment No Do not ingest | No data avaii IA No data xposure poprtively ious amounts Handling a for Spilled water. | Skin? Yes lable for othe RC Monographic s of water. and Use respirator with skin. even with skin. even | Yes r routes phs? C No No Special Other Special Other StiorSplash process and clothing | Ingestion? DSHA Regulation? data DSHA Regulation? data or cartridge. Yes None of goggles Wash thoroughly |

| EDVOTEK. | May be Standard | Mat used to I. 29 CF | comply with OSHA's R 1910.1200 Standa specific requirement | ata Sheet Hazard Communi ard must be consu ents. | ication Ilted for |
|--|--------------------|----------------------------|--|---|----------------------|
| IDENTITY (As Used on Label and List) | | | Note: Blank spaces a | re not permitted. If | any item is not |
| InstaStain® Ethidium Bromid | e | | be marked to indicate | that. | the space must |
| Section I | | | | | |
| Manufacturer's Name | | Emer | gency Telephone N | umber (301) 2 | 251-5990 |
| InstaStain, Inc. | | Teleni | one Number for infor | mation | |
| P.O. Box 1232 | | 10.00 | | (301) 2 | 51-5990 |
| West Bethesda, MD 20 | 827 | Date | Prepared | 11/21/10 |) |
| | | Signa | ure of Preparer (optic | nal) | |
| Section II - Hazardous Ingred | lients/Iden | tify Ir | formation | | |
| Hazardous Components [Specific Chemical Identity; Common Name(s)] | OSHA | PEL | ACGIH TLV R | Other Limits ecommended | % (Optional) |
| Ethidium Bromide | 1 | Data no | t available | | |
| (2,7-Diamino-10-Ethyl-9-Phenyl | phenanthridin | ium Bro | omide) | | |
| CAS# 139-33-3 | | | | | |
| Section III - Physical/Chemic | al Charact | eristi | cs | | |
| Boiling Point | No data | Spe | cific Gravity ($H_20 =$ | 1) | No data |
| Vapor Pressure (mm Hg.) | No data | Mel | ting Point | | No data |
| Vapor Density (AIR = 1) | No data | Eva (Bu | poration Rate tyl Acetate = 1) | | No data |
| Solubility in Water Soluble | | | | | |
| Appearance and Odor Chemical bour | nd to paper, no | o odor | | | |
| Section IV - Physical/Chemic | al Charac | teristi | cs N.D. = | No data | |
| Flash Point (Method Used) No d | ata | Flam | mable Limits | LEL N.D. | UEL N.D. |
| Extinguishing Media Water spray, ca | arbon dioxide, | dry ch | emical powder, alco | hol or polymer i | foam |
| Special Fire Fighting Procedures | Vear protective | e clothi | ng and SCBA to pre | vent contact wit | h skin & eyes |
| Unusual Fire and Explosion Hazards | Emits toxic fur | nes | | | |

| Section V - Reactiv | | | | | | | |
|--|---|-------------------------|---|--------------------|----------------------|--|--|
| Stability | Unstable | Unstable Condit | | | itions to Avoid | | |
| | Stable | Х | None | | | | |
| Incompatibility | Strong oxidi | zing agents | 3 | | | | |
| Hazardous Decomposition | or Byproducts | | | | | | |
| Carbon mo | onoxide, Carbon did | xide, nitroc | en oxides, hydrog | en bromide | e gas | | |
| Hazardous Polymerization | May Occur | | Conditions to Avoid | | | | |
| | Will Not Occur | Х | None | | | | |
| Section VI - Healti | h Hazard Data | | | | | | |
| Houte(s) of Entry: | Inhalatio | on? Yes | Skin? | Yes | Ingestion? Yes | | |
| Health Hazards (Acute a Acute: Material i | and Chronic) Chror rritating to mucous i | nic: May al membrane | er genetic materia s, upper respirator | ıl y tract, eye | s, skin | | |
| Carcinogenicity: No dat | a available NTP? | | IARC Monogra | phs? | OSHA Regulation? | | |
| Signs and Symptoms of | Exposure Irritati | on to muce | us membranes and | l upper res | piratory tract | | |
| Medical Conditions Gen | erally Aggravated b | y Exposure | No data | | | | |
| | | | | | | | |
| Emergency First Aid Pro | Dcedures | | i solles and assessed | | | | |
| | Treat | symptoma | ically and support | Ivery | | | |
| - | | | | | | | |
| Section VII - Preca | utions for Safe | e Handli | ng and Use | | | | |
| Steps to be Taken in cas | e Material is Releas | sed for Spil | led | | | | |
| Wear SC | BA, rubber boots, r | ubber glov | es | | | | |
| Waste Disposal Method | Mix material wi | ith combust | ible solvent and b | urn in a che | emical incinerator | | |
| | equipped afterb | urner and s | crubber | | | | |
| | | | | | | | |
| Precautions to be Taken | in Handling and St | oring | | | | | |
| | Use in chemical | fume hood | with proper prote | ective lab g | ear. | | |
| Other Precautions | | | | | | | |
| Other Trecadions | Mutagen | | | | | | |
| | | | | | | | |
| Section VIII - Cont | rol Measures | | | | | | |
| Respiratory Protection (| Specify Type) | SCBA | | | | | |
| Ventilation | Local Exhaust | Yes | | Special | Chem. fume hood | | |
| | Mechanical (Ge | eneral) ¹ | No | Other | None | | |
| Protective Gloves R | ubber | | Eye Protec | tion | Chem. safety goggles | | |
| Other Protective Clothing | g or Equipment | Rubber boo | ots | | | | |
| Work/Ukraina Drastian | | | | | | | |
| work/Hygienic Practices | 5 | Use in cher | nical fume hood w | ith proper | protective lab gear. | | |

| ED VOTEK. | May be used Standard. 29 | Interial Safety Data Sheet d to comply with OSHA's Hazard Communication CFR 1910.1200 Standard must be consulted for specific requirements. | ı or | | |
|--|-----------------------------|--|----------------------------------|--|--|
| IDENTITY (As Used on Label and List) | | Note: Blank spaces are not permitted. If applicable, or no information is available, | any item is not the space mus | | |
| InstaStain® Blue, Methylene Blue I | Plus™ | be marked to indicate that. | | | |
| Section I | | Farmer Televisor Number | | | |
| Manufacturer's Name | | (301) 251-5990 | | | |
| EDVOTEK, Inc. | | Telephone Number for information | | | |
| Address (Number, Street, City, State, | Zip Code) | (301) 2 | 51-5990 | | |
| 14676 Rothgeb Drive | | Date Prepared 11-21-10 | | | |
| | | Signature of Preparer (optional) | | | |
| Section II - Hazardous Ingred | lients/Iden | ntify Information | | | |
| Hazardous Components [Specific Chemical Identity; Common Name(s)] | OSHA | Other Limits A PEL ACGIH TLV Recommended | % (Optiona | | |
| Methylene Blue | | | | | |
| 3.7 Bis (Dimethylamino) Phenothiazin | 5 IUM Chlo | oride No data available | | | |
| CAS # 61-73-4 | | | | | |
| Section III - Physical/Chemic | al Charact | teristics | | | |
| Boiling Point | No data | Specific Gravity ($H_2^0 = 1$) | No data | | |
| Vapor Pressure (mm Hg.) | No data | Melting Point | No data | | |
| Vapor Density (AIR = 1) | No data | Evaporation Rate (Butyl Acetate = 1) | No data | | |
| Solubility in Water Soluble - cold | | | | | |
| Appearance and Odor Chemical bou | nd to paper, n | 10 odor | | | |
| Section IV - Physical/Chemic | al Charac | teristics | | | |
| Flash Point (Method Used) No data ava | ulable | Flammable Limits LEL | UEL No data | | |
| Extinguishing Media Water spray | , carbon dioxi | ide, dry chemical powder, alcohol or polymo | er foam | | |
| Special Fire Fighting Procedures Self contained breathing app | aratus and pro | otective clothing to prevent contact with skin | n and eyes | | |
| Unusual Fire and Explosion Hazards | Emits toxid | fumes under fire conditions | | | |

| Stability | Unstable | | Conditions to Avoid | |
|---|--------------------------------|-----------|------------------------------|-----------------------------|
| | Stable | Х | None | |
| Incompatibility Stro | ng oxidizing agen | its | | |
| Hazardous Decomposition or I Toxic fumes of Carbon n | Byproducts nonoxide, Carbon | dioxide | nitrogen oxides, sulfur oxi | des, hydrogen, chloride gas |
| Hazardous | May Occur | | Conditions to Avoid | , |
| Polymerization | Will Not Occur | X | None | |
| Section VI - Health H | azard Data | | | |
| Route(s) of Entry: | Inhalatio | n? Yes | Skin? Yes | Ingestion? Yes |
| Health Hazards (Acute and Skin: May cause ski | Chronic) n irritation Ev | es: May | v cause eve irritation In | halation: Cyanosis |
| Carcinogenicity: | NTP? | | IARC Monographs? | OSHA Regulation? |
| Meets criteria for pro | posed OSHA med | dical rec | ords rule PEREAC 47.3042 | 0.82 |
| Signs and Symptoms of Ex | posure No da | ta availa | ble | |
| Medical Conditions Genera | lly Aggravated by | Exposu | re No data available | |
| Emergency First Aid Proce | dures Treat syr | nptomat | ically | |
| Section VII - Precaut | ions for Safe | Hand | ling and Use | |
| Steps to be Taken in case N | laterial is Releas | ed for Sp | pilled | |
| | Ventilate | area an | d wash spill site | |
| Waste Disposal Method | Mix mate | erial wit | h a combustible solvent and | l burn in chemical |
| incinerator equipped wit | h afterburner and | scrubbe | r. Check local and state reg | ulations. |
| Precautions to be Taken in | Handling and Sto | ring | | |
| Keep tightly closed. Sto | re in cool, dry pla | ice | | |
| Other Precautions | | | | |
| None | | | | |
| Section VIII - Contro | I Measures | | | |
| Respiratory Protection (Sp | ecify Type) MI | IOSH/O | SHA approved, SCBA | |
| Ventilation | Local Exhaust | | Special | |
| | Mechanical (Ge | neral) | Required Other | |
| Protective Gloves Rubb | er | | Eye Protection | Chem. safety goggles |
| | | | | |
| Other Protective Clothing o | r Equipment | Rubber | boots | |