

EDVOTEK® Quick Guide: Agarose Gel Electrophoresis



© All rights reserved, Edvotek, Inc. 2018

What is Electrophoresis?

Electrophoresis is a technique that allows us to separate DNA, RNA or proteins according to their size.

What do I need to separate a mixture of DNA molecules?

In addition to your DNA sample, you will need:

- Gel Loading Solution – includes glycerol to help DNA samples enter into the wells and a visible dye to monitor migration through the gel.
- Agarose – a polysaccharide used as the separation matrix.
- Electrophoresis Buffer – contains ions necessary to conduct an electrical current, maintains pH of experiment.
- Horizontal electrophoresis apparatus – holds the buffer and the gel, has positive and negative electrodes.
- Power supply – generates the current necessary to move DNA through gel.
- Micropipet – used to transfer samples into wells.
- A special stain that allows us to visualize DNA.

How does electrophoresis separate DNA fragments?

The mixture of DNA molecules is added into depressions (or “wells”) within a gel, and then an electrical current is passed through the gel (Figure 1A). Because the sugar-phosphate backbone of DNA has a strong negative charge, the current drives the DNA through the gel towards the positive electrode (Figure 1B).

At first glance, an agarose gel appears to be a solid at room temperature. On the molecular level, the gel contains small channels through which the DNA can pass. Small DNA fragments move through these holes easily, but large DNA fragments have a more difficult time squeezing through the tunnels. Because molecules with dissimilar sizes travel at different speeds, they become separated and form discrete “bands” within the gel. After the current is stopped, the bands can be visualized using a stain that sticks to DNA (Figure 1C).








Figure 1: Overview of Agarose Gel Electrophoresis



EDVOTEK® Quick Guide: Agarose Gel Electrophoresis







1. **DILUTE** concentrated (50X) buffer with distilled water to create 1X buffer (see Table A).
2. **MIX** agarose powder with 1X buffer in a 250 ml flask (see Table A).
3. **DISSOLVE** agarose powder by boiling the solution. **MICROWAVE** the solution on high for 1 minute. Carefully **REMOVE** the flask from the microwave and **MIX** by swirling the flask. Continue to **HEAT** the solution in 15-second bursts until the agarose is completely dissolved (the solution should be clear like water).
4. **COOL** agarose to 60° C with careful swirling to promote even dissipation of heat.
5. While agarose is cooling, **SEAL** the ends of the gel-casting tray with the rubber end caps. **PLACE** the well template (comb) in the appropriate notch.
6. **POUR** the cooled agarose solution into the prepared gel-casting tray. The gel should thoroughly solidify within 20 minutes. The gel will stiffen and become less transparent as it solidifies.
7. **REMOVE** end caps and comb. Take particular care when removing the comb to prevent damage to the wells.

Size of Gel Casting tray	Concentrated Buffer (50x)	+ Distilled Water	+ Amt of Agarose	= TOTAL Volume
7 x 7 cm	0.6 ml	29.4 ml	0.23 g	30 ml
7 x 10 cm	1.0 ml	49.0 ml	0.39 g	50 ml
7 x 14 cm	1.2 ml	58.8 ml	0.46 g	60 ml



© All rights reserved, Edvotek, Inc. 2018

EDVOTEK® Quick Guide: Agarose Gel Electrophoresis



8. **PLACE** gel (on the tray) into electrophoresis chamber. Completely **COVER** the gel with 1X electrophoresis buffer (See Table B for recommended volumes).
9. **LOAD** entire sample volumes into wells in consecutive order.
10. **PLACE** safety cover. **CHECK** that the gel is properly oriented. Remember, the samples will migrate toward the positive (red) electrode.
11. **CONNECT** leads to the power source and **PERFORM** electrophoresis (See Table C for time and voltage guidelines).
12. After electrophoresis is complete, **REMOVE** the gel and casting tray from the electrophoresis chamber and proceed to **STAINING & VISUALIZATION**.

Table B

1x Electrophoresis Buffer (Chamber Buffer)

EDVOTEK Model #	Total Volume Required	50x Conc. Buffer	Dilution + Distilled Water
M6+ & M12 (new)	300 ml	6 ml	294 ml
M12 (classic)	400 ml	8 ml	392 ml
M36	1000 ml	20 ml	980 ml



Wear gloves and safety goggles

Table C

Time and Voltage Guidelines (0.8% Agarose Gel)

Volts	Electrophoresis Model		
	M6+	M12 (new)	M12 (classic) & M36
	Min. / Max.	Min. / Max.	Min. / Max.
150	15 / 20 min.	20 / 30 min.	25 / 35 min.
125	20 / 30 min.	30 / 35 min.	35 / 45 min.
75	35 / 45 min.	55 / 70 min.	60 / 90 min.

Reminder:

Before loading the samples, make sure the gel is properly oriented in the apparatus chamber.

Electrophoresis Equipment

See the **EQUIPMENT** section in our Resource Guide for our full range of electrophoresis and power supplies or visit our website at: www.edvotek.com



Cat. #515
**M36 HexaGel™ DNA
Electrophoresis Apparatus**
For 6 Lab Groups



Cat. #502/504
**M12 Complete™
Electrophoresis Package**
For 2 Lab Groups



Cat. #500
**M6Plus Electrophoresis
Apparatus**
For 1 Lab Group



Cat. #5010
TetraSource™ 300
30-300 V for 1 to 4 units



Cat. #509
DuoSource™ 150
75/150 V for 1 or 2 units



Cat. #507
DuoSource™ 75
75 V for 1 or 2 units



Cat. #589-#593
EDVOTEK® Variable Micropipets
From 0.1 µl to 5000 µl



Cat. #585-#588
EDVOTEK® Fixed Volume MiniPipets™
From 5 µl to 200 µl



Cat. #541
EdvoCycler™
Holds 25 x 0.2 ml tubes



Cat. #558
Midrange UV Transilluminator
7x14 cm UV filter.