

CARE AND FEEDING OF *C. ELEGANS*

| PROBLEM: | CAUSE: | ANSWER: |
|---|---|--|
| Received source plates and worms do not move when observed through a microscope | Worms could die from extreme temperatures during shipping. | Acquire new plate of worms. |
| | Worms are dormant. | Feed the worms by adding one BactoBead™ to 1ml of LB or water. Grow the bacteria for one hour at 37°C. Add 200 µl of the culture to each plate. |
| Stored the source plates at 4°C for a week, and they don't move. | Worms are dormant. | Feed the worms by adding one BactoBead™ to 1ml of LB or water. Grow the bacteria for one hour at 37°C. Add 200 µl of the culture to each plate. |
| While transferring the worms to a new plate, the agar was flipped in the wrong orientation. | It is difficult to transfer a small piece of agar. | The direction of the agar does not matter, worms will leave the chunk and migrate to their food. |
| Can not see the worms under the microscope after transfer. | Worms require a short amount of time to migrate to plate. | Wait for 5 - 10 minutes and observe under the microscope again. |
| | Chunk did not contain enough live worms. | Make new chunk of agar from source plate and repeat the transfer process onto the plates. |
| The chunked plates containing the worms show a bright white or pink growth. | The plates are contaminated. | Prepare a new plate of worms from a clean source plate. |
| | | Transfer a chunk of agar from a clean region of the plate to a fresh NGM agar plate. |
| Do not see growth of OP50 on plate after incubation. | NGM is a clear agar and OP50 bacteria forms a limited lawn. | Using a loop gently scrape the top of the agar to check for bacterial growth or smell plate for bacterial odor. If no growth is observed, plate more OP50 bacteria. |
| Can not count the worms under the microscope. | The microscope is set to higher magnification. | Start with the lower magnification and try to visualize many locations around the plate. |
| | Too few worms on plate. | Wait three more days before performing the experiment. During this time, be sure to feed the worms every few days by adding one BactoBead™ to 1ml of LB or water. Grow the bacteria for one hour at 37°C. Add 200 µl of the culture to each plate. |
| Can not see worms in the counting chamber. | Agar from plate has contaminated worms during washing. | Avoid removing agar from the plate by gently pipetting and avoiding contact with the surface. |
| | Worms are at the bottom of the tube. | Gently flick the tube to resuspend worms prior to pipetting. |
| Need to dispose of the worms and NGM plates. | Experiment has been completed. | Plates should be soaked in 10% bleach and discarded according to school regulations. |
| After chunking the worms and left them at RT, the NGM agar is too thin or does not cover the bottom of the plate. | Incorrect volume was poured into the plates. | It is important to have the correct thickness of agar in the plates. Be sure to pour the amount recommended by the protocol. |
| | The plates were left uncovered and the agar dehydrated. | Keep plates inside of the box. If it is very hot in your classroom, seal plates with parafilm. |
| Received worms two or more weeks before performing the experiment. | Worms may be dehydrated or starving. | Pour the NGM plates, growth bacteria and chunk the plates into a new plate as soon as possible. |