

EDVOTEK® INQUIRY GUIDE:
ELECTROPHORESIS



TEACHERS:

For the answer key, please
contact us at info@edvotek.com

www.edvotek.com

EDVOTEK®

Key Terms

DEFINE the following terms before you begin the guide then refer to them as needed. Use any resources you want!

Electrophoresis:

Matrix:

Agarose:

DNA:

RNA:

Protein:

Electrical Field:

Gel:


PCR:

Restriction Enzyme:

Ladder:



- This question requires a pencil and a thicker marker. **TRAVERSE** the matrix below with your pencil, moving exclusively through the white areas and steering clear of the black dots. Start from the left and proceeding to the right. **REPEAT** with your marker. Finally, **ANSWER** the four question on the right.




- Which run was faster?
- How does this relate to electrophoresis?
- What does the pencil represent in this simulation?
- What does the marker represent in this simulation?


- CIRCLE** the molecules below that would work on the electrophoresis apparatus below.

Molecule:	DNA (- charge) (Lane 1)	RNA (- charge) (Lane 2)	Protein (neutral) (Lane 3)	Protein (+ charge) (Lane 4)	Protein (- charge) (Lane 5)
------------------	----------------------------	----------------------------	-------------------------------	--------------------------------	--------------------------------


- In the illustration below, **DRAW** an arrow (<← or →) in each molecule's corresponding "lane" to indicate which direction (if any) each molecule will move.



4. In DNA electrophoresis, the matrix used to contain and separate molecules is called a gel. Typically, this gel is made from agarose. By adding more or less agarose, scientists can separate a wide range of DNA molecules from the very big (25,000 bp) to the very small (<10 bp).
- a. **CALCULATE** how much agarose you'd need to create 50 mL of a 1% gel using the weight/volume formula below. (A 1% gel is optimal for separating DNA fragments 500-10,000 bp long.)
- Agarose (g) = Percentage (expressed as a decimal or fraction) x Gel Volume (mL)
- b. If you were working with DNA fragments between 50-500 bp, would you make a lower or higher percentage gel (i.e. add less or more agarose)?
5. Congratulations! You've successfully PCR amplified your DNA samples and run the results on a gel. Now it's time to analyze the results. Use the gel below to **ANSWER** the following questions. (You will need access to colored pens.)
- a. In RED circle the two largest DNA fragments in the three samples.
- b. In BLUE circle the two smallest DNA fragments in the three samples.
- c. In BLACK circle DNA fragments larger than 500 base pairs in the three samples.



6. Congratulations again! You've successfully amplified your DNA samples, completed a restriction enzyme reaction*, and run the results on a gel. Now it's time to analyze the results. Use the gel below to **COMPLETE** the table.



Lane	Number of Fragments	Estimated Size of Fragments	Estimated Sum of all Fragments
Marker	7	6751, 3652, 2827, 1568, 1118, 825, 630	---
Sample 1 (Undigested DNA)	1	7720	7720
Sample 2 (DNA digested with EcoRI)			7720
Sample 3 (DNA digested with BamHI)			7720
Sample 4 (DNA digested with EcoRI and BamHI)			7720


* During a restriction enzyme reaction amplified DNA is mixed with enzymes that cut the DNA at specific DNA sequences resulting in smaller fragments. The more frequently a specific sequence appears the more often the DNA will be cut and the more fragments.

7. Now do the reverse. Use the blank agarose canvas and table below to **DRAW** the results from a diagnosis experiment. The marker in Lane 1 is already filled in for you.

Lane	Number of Fragments	Estimated Size of Fragments
Marker (already filled in)	7	6751, 3652, 2827, 1568, 1118, 825, 630
Negative Control	1	4282
Positive Control	3	4282, 3000, 1282
Patient 1	3	4282, 3000, 1282
Patient 2	1	4282



8. Use this next agarose canvas to **DRAW** a simple shape out of bands (e.g. a circle, square, heart, emoji). What base pair bands would you need to make this shape in the lab?



Lane	Number of Fragments	Estimated Size of Fragments
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
11		
12		

NOTE: Remember that in electrophoresis the key measurement is distance traveled from wells so feel free to use different colors, thicknesses, or textures for the bands.