

EDVOTEK

PCR Perfection



PLTW SUMMIT



www.EDVOTEK.com

 Follow @Edvotek



EDVOTEK

The **Biotechnology Education** Company




Celebrating 30 years of science education!

EDVOTEK


The **Biotechnology Education** Company



Experiments



Reagents



Equipment



Resources

Edvotek PLTW Experiments

Kit	Description
118-PLTW → 418	Hypercholesterolemia Kit
225-PLTW → 425	DNA Detectives Kit
235-PLTW → 435	DNA Microarrays
268-PLTW → 468	Mystery Infections Kit
303 - PLTW → 403	Exploring Biotechnology with GFP
339 - PLTW → 439	Sequencing The Human Genome Kit
990-PLTW → 490	Morphology of Cancer Cells Kit
953	Multiplex PCR Testing of Water Contaminants
345	SNP Analysis of the PTC Gene
491 (1.1.4)	Forensic Blood Typing
401 (4.2.3)	Under the Sea
485	Preparation of Metaphase Chromosome Spread
430	DNA Fingerprinting



Edvotek PLTW Equipment

Kit	Description
504	M12 Complete™ Electrophoresis Apparatus
5010	TetraSource™ 300 Digital Electrophoresis PS
557	TruBlu™ Blue Light Transilluminator
589 / 591 / 592-1	Edvotek Variable Micropipets, 0.5-10 µl, 2-20 µl, 10-100 µl, 100-1000 µl
540 / 541-542	Edvocycler Jr, Edvocycler 2



Today's Strategy

New water quality extraction - Tips & Tricks for PTC



Water Quality Testing Extraction

- Updated the extraction protocol this summer
- Removed incubation steps
- Check literature link on freezer bag
- Will still work with old components



Water Quality Testing Extraction


Preparation For:	What to do:	When:	Time Required:
Module I: Isolation of Bacterial DNA from Contaminated Water	Prepare contaminated water samples and DNA extraction solutions.	No more than one hour before class.	45 min.
	Prepare and aliquot Lysis Buffer.	Prepare on the day the students will be performing the experiment OR freeze for up to one week.	15 min.
	Preheat water baths.	Anytime before performing the experiment.	15 min.
Module II: Multiplex PCR Amplification of Water Contaminants	Program thermal cycler.	Anytime before performing the experiment.	15 min.
	Prepare and aliquot various reagents (Primer, DNA template, ladder, etc.)	One day to 30 minutes before performing the experiment.	30 min.
Module III: Separation of PCR Product by Electrophoresis	Prepare 1X Electrophoresis buffer and dilute SYBR® Safe.	Up to one week before performing the experiment.	45 min.
	Prepare molten agarose and pour gel.		
Module IV: Staining Agarose Gels with FlashBlue™ (OPTIONAL)	Prepare staining components.	The class period or overnight after the class period.	10 min.

Isolation of DNA from Bacteria

FOR MODULE I


Each Group should receive:

- 4 0.5 mL screw-cap tubes containing contaminated water samples (one of each)
- 4 1.5 mL screw-top tubes containing Lysis Buffer, on ice



- Step 1: Label tubes with bacteria
- Step 2: Add 100 μ L water/bacteria sample to lysis buffer
- Step 3: Incubate at 55 deg. C for 5 min

Isolation of DNA from Bacteria



- Step 5: Vortex or flick to mix
- Step 6: Incubate at 99 deg. C for 5 min
- Step 7: Place on ice and proceed to PCR

Have the old freezer bag? No problem!

Components	Storage	Check (✓)
• PCR EdvoBeads™ <i>Each PCR EdvoBead™ contains dNTP mixture, Taq DNA polymerase buffer, Taq DNA polymerase, MgCl₂, and reaction buffer</i>	Room Temp., desiccated	<input type="checkbox"/>
A Universal DNA Buffer	-20° C Freezer	<input type="checkbox"/>
B TE buffer	-20° C Freezer	<input type="checkbox"/>
C LyphoPrimer™ Mix <i>Contains E. coli-specific primers, Bacillus subtilis-specific primers, and Serratia marcescens-specific primers</i>	-20° C Freezer	<input type="checkbox"/>
D LyphoControl™ (Complete PCR Control)	-20° C Freezer	<input type="checkbox"/>
E EdvoQuick™ DNA ladder	-20° C Freezer	<input type="checkbox"/>
F Proteinase K	-20° C Freezer	<input type="checkbox"/>
G Potassium Acetate	-20° C Freezer	<input type="checkbox"/>
H DNA Extraction Buffer	-20° C Freezer	<input type="checkbox"/>
• E. coli BactoBeads™	4° C, desiccated	<input type="checkbox"/>
• Bacillus subtilis BactoBeads™	4° C, desiccated	<input type="checkbox"/>
• Serratia marcescens BactoBeads™	4° C, desiccated	<input type="checkbox"/>

This experiment is designed for 6 groups of students.

LyphoPrimer™

LyphoControl™

IMPORTANT NOTE:

This protocol has been optimized to reduce time and improve results. Because of this, *components G and H will no longer be used.*

Amplification of the bacteria

Teacher preparation:

Preparation For:	What to do:	When:	Time Required:
Module II: Amplification of the PTC Region	Prepare and aliquot various reagents (Primer, DNA template, ladder, etc.)	One day to 30 min. before performing the experiment.	30 min.
	Program Thermal Cycler	One hour before performing the experiment.	15 min.

- Student experiment should take 15-20 minutes.
- PCR reactions should be setup and run on same day.
- PCR products can be frozen and stored until needed.

DNA Amplification using PCR

- DNA template
- Primers
 - Short piece of DNA that defines the area to amplify
- Thermostable DNA polymerase (*Taq*)
- Free nucleotides (dNTPs)
- Buffer




DNA Amplification using PCR

- DNA template
- Primers
 - Short piece of DNA that defines the area to amplify
- Thermostable DNA polymerase (*Taq*)
- Free nucleotides (dNTPs)
- Buffer




DNA Primers


- 50 kb segment of Lambda Phage DNA



PCR Amplifies Specific DNA Sequences




PCR Amplifies Specific DNA Sequences



- **Denaturation (94°C)**
Double-stranded DNA is “unzipped” into single strands.
- **Annealing (40-65°C)**
Primers base pair with the target DNA sequence.
- **Extension (72°C)**
Taq polymerase extends the primer and synthesizes a new strand of DNA.


PCR Amplifies DNA Exponentially



PCR Amplifies DNA Exponentially



Electrophoresis at a Glance




Electrophoresis Separates DNA Fragments By Size

- The sugar-phosphate backbone of DNA has a strong negative charge.
- When an electrical current is passed through the gel, the current drives the DNA fragments through the gel towards the positive electrode.




Electrophoresis Separates DNA Fragments By Size

- The gel contains small channels through which the DNA can pass.
- Small DNA fragments move through these holes easily, but large DNA fragments have a more difficult time squeezing through the tunnels.



Electrophoresis Separates DNA Fragments By Size

- Because molecules of different sizes travel at different speeds, discrete bands are formed.
- After the current is stopped, the bands can be visualized using a stain that sticks to DNA.
 - InstaStain® EtBr
 - InstaStain® Blue
 - FlashBlue™ Stain
 - Sybr® Safe stain




SYBR® Safe DNA Stain

In-gel Staining

1. Melt agarose and cool to 65°C.
2. Add concentrated Sybr® Safe stain to the molten gel at 1:10,000 dilution (5 µL per 50 mL agarose solution).
3. Run DNA samples through gel – no post staining or destaining necessary!

Post-electrophoresis Staining

1. Dilute concentrated stain to 1:20,000 (5 µL per 100 mL distilled water).
2. After electrophoresis, place gel in tray. Cover gel with diluted Sybr® Safe stain.
3. Stain gel for 10 – 15 minutes.



Kit #109
Transilluminator #558
SybrSafe® Stain #608

TruBlu2™ BlueLight Transilluminator

- Optimized for **SYBR® Safe** stained gels • Option for blue or white light •
- Large viewing area • No harmful UV •



1.800.EDVOTEK
www.edvotek.com



Multiplex PCR Testing of Water Contaminants



Lane	Sample
1	EdvoQuick™ DNA Ladder
2	Extracted DNA - E. coli
3	Extracted DNA - Serratia
4	Extracted DNA - B. subtilis
5	Extracted DNA - Serratia & E. coli
6	Control DNA - E. coli

EDVO-Kit #953

REVISED
&
UPDATED

Edvo-Kit #

345

Edvo-Kit #345

Exploring the Genetics of Taste: SNP Analysis of the PTC Gene Using PCR

Experiment Objective:

In this experiment, students will isolate their own DNA and use PCR to amplify a segment of the *TAS2R38* gene, which is responsible for detecting the bitter taste of PTC. Digestion of the PCR products and analysis by agarose gel electrophoresis are used to differentiate tasters and non-tasters. Finally, students' genotype is linked to phenotype by tasting the PTC paper.

See page 3 for storage instructions.

LyphaPrimer™

LyphaControl™



PTC Improvements

- Lyophilized reagents increase stability
- Multiple stopping points
- Complete Control
- Larger volumes
- PCR cycling time reduced
- Instructor guide

Exploring the Genetics of Taste: SNP Analysis of the PTC Gene Using PCR EDV-Kit 345

Strategy & Planning

Before you start:

- Ensure that the lot version (S45, 192514) matches the label found in your peripherals.
- Review the experiment background with students.
- Remind students of safety guidelines and disposal information for their cheek cell solutions.

Techniques:
Before beginning this experiment, students should be proficient in pipetting.

Pre-Lab Preparation

Detailed pre-lab preparation instructions can be found on page 20. Remember to plan for stopping points if you will not be continuing directly to Module II.

IMPORTANT:
Only use saline for the student cheek cell rinse. Using sports drinks or other solutions may inhibit the PCR in Module II.

Resources

For additional information about this experiment, visit www.edvotek.com/S45

For resources on pipetting, isolation of cheek cell DNA, and more: www.edvotek.com/Quick-Guides and www.youtube.com/EdvotekInc

STUDENT PROTOCOL PREVIEW

Module I Overview

In Module I, you will isolate DNA from your cheek cells. First, you will vigorously rinse your mouth with saline (salt water), which will dislodge cells into the solution. The cells are gathered using a centrifuge to pellet them at the bottom of a microcentrifuge tube, allowing the saline to be removed. Next, a lysis buffer is added and the solution is incubated at 55° C and 95° C to break open the cells and release the DNA. Finally, the cell lysis is centrifuged: this will collect the cell debris in a pellet while leaving the DNA in the supernatant. The DNA-containing supernatant will be used in Module II.

Experiment Overview

Where's my DNA?

1. Collecting cheek cells (5 min)

2. Cell lysis (5 min)

3. Isolated DNA (5 min)

4. After incubation, Step 4 and 5: Incise the cell pellet

5. After lysis, Step 6: Mix in the lysis buffer

6. After centrifugation, Step 7: In the supernatant

7. Cell debris

8. Isolated DNA

EDVOTEK, Inc. 1.800.EDVOTEK • Fax: 202.703.5351 • info@edvotek.com • www.edvotek.com

10

In-Depth Understanding of Module I:

In Module I, the students will be isolating DNA from cheek (Buccal) cells using a quick cell lysis protocol. The cells are collected in a saline solution that has been prepared in advance - it is essential that saline is used here, since other solutions can inhibit the PCR in Module II. Once collected, the cheek cells are suspended in a lysis buffer containing a pH buffer, detergent, and protease. This solution, as well as two brief incubations, will weaken and then lyse the cells. Finally, a DNA solution is collected by centrifuging the cell lysate. It is important that students collect a relatively pure DNA solution as cell debris from the lysis steps can inhibit PCR. It might be helpful to review the "Where's my DNA" flowchart before beginning the experiment to ensure that students do not accidentally discard their DNA.

10 Duplications of any part of this document is permitted for non-profit educational purposes only. Copyright © 2015-2019 Edvotek, Inc., all rights reserved. S45-19051-E-1



Instructor's Guide

STUDENT PROTOCOL PREVIEW

Module I: Isolation of DNA from Human Cheek Cells



- LABEL** an empty 1.5 mL screw top microcentrifuge tube and a cap of saline with your lab group and/or initials.
- RINSE** your mouth vigorously for 30 seconds using 30 mL saline solution. **EXPEL** the solution back into the same cup.
- SWIRL** the cup gently to resuspend the cells. **TRANSFER** 1.5 mL of the cell solution into the tube with your initials.
- CENTRIFUGE** the cell suspension for 2 minutes at full speed to pellet the cells. **POUR** off the supernatant (the liquid above the cell pellet) (do **NOT** disturb the cell pellet).
- RESUSPEND** steps 3 and 4 once more.
- RESUSPEND** the cheek cell pellet in 150 µL lysis buffer by pipetting up and down or by vortexing vigorously. **NOTE:** ensure that the cell pellet is fully resuspended and that no clumps of cells remain.
- CAP** the tube and **PLACE** it in a water bath that **INCUBATE** the sample in a 55°C water bath for 5 minutes.
- RINSE** the sample by vortexing or by flicking the tube vigorously for 20 seconds.
- INCUBATE** the sample in a 99°C water bath for 5 minutes. **NOTE: STUDENTS MUST USE SCOVING TIPS WHEN CHANGING SAMPLES.**
- CENTRIFUGE** the cell lysate for 2 minutes at full speed.
- TRANSFER** 50 µL of the supernatant to a clean, labeled microcentrifuge tube. **PLACE** the tube in ice.

The extracted DNA is now ready for Module II: Amplification of the PTC Region. If you are ready to proceed, turn to page 12. Alternatively, the extracted DNA may be stored in the **FREEZER** (-20°C) until needed.

OPTIONAL STOPPING POINT:
The extracted DNA may be stored in the freezer (-20°C) until needed.

1.800.EDVOTEK • Fax: 201.770.1581 • info@edvotek.com • www.edvotek.com

Small block of any part of this document is permitted for non-profit educational purposes only. Copyright © 2015-2019 EDVOTEK, Inc., all rights reserved. 345-1905(1.1)

EDVOTEK

11

Vocabulary Review

Saline - A salt water solution.
Resuspend - Place a solid back into solution.
Supernatant - The liquid layer above a sediment or solid.
Incubate - To maintain a solution at an ideal condition or temperature.

Experimental Assessment

Ask the students to critically analyze each step of the protocol. For example: Why is the saline rinse centrifuged twice, or why do we only collect 50 µL of the supernatant in step 11?

Stopping Point

DNA is very stable when frozen. Once the students have completed this module, their samples can be stored in the freezer (-20°C) until needed in Module II.

Safety

Human saliva may transmit disease. We recommend sterilizing all human samples with 10% bleach after completing the experiment.

Remind students to wear gloves and goggles while working in the lab, and to leave the lab area before eating or drinking (after the saline rinse).

Boiling water is hot, use caution!


Strategy & Planning

- You must use saline for the rinse!
- After centrifuging you should see at least a faint pellet in each student's tube. If not, ask the students to swirl their remaining saline solution to resuspend settled cells and repeat the steps 3 and 4 again.
- Ensure cells are fully resuspended - no clumps!
- This does not need to be a rolling boil, but we want the water as hot as possible.
- Ensure that there is no debris transferred with the DNA supernatant. If you notice any clumps or cloudiness in the samples they can be centrifuged again to pellet any debris.

STUDENT PROTOCOL PREVIEW

Module II: Amplification of the PTC Region

Now that you have isolated your DNA, during Module I, the next step is to amplify a specific region of the *UGT2C8* gene. First, you will combine your DNA (red) with a mixture of PCR primers (yellow) and a PCR Inhibitor™ PLUS, creating a PCR sample. Once the sample has been prepared it will be placed into a thermal cycler and the DNA will be amplified by PCR.



- COMBINE** the red extracted DNA from Module I.
- LABEL** a fresh 0.2 mL PCR tube with your initials.
- ADD** 40 µL PCR primer mix (yellow), 10 µL extracted DNA (red), and a PCR Inhibitor™ PLUS.
- MIX** the PCR sample. Shake using the PCR Inhibitor™ PLUS is completely dissolved. If mixed correctly, the final solution will be light orange.
- CENTRIFUGE** the tube for a few seconds to collect the sample at the bottom of the tube.
- AMPLIFY** the DNA using PCR. **PCR cycling conditions:**
 - Initial denaturation 94°C for 4 minutes
 - 94°C for 30 seconds
 - 45°C for 30 seconds
 - 72°C for 30 seconds
 - Final extension 72°C for 5 minutes

The amplified DNA is now ready for Module III: Restriction Digest of the PCR PCR product. If you are ready to proceed, turn to page 13. Alternatively, the amplified DNA may be stored in the **FREEZER** (-20°C) until needed.

OPTIONAL STOPPING POINT:
The PCR samples may be stored at -20°C for restriction digest at a later time.

1.800.EDVOTEK • Fax: 201.770.1581 • info@edvotek.com • www.edvotek.com

Small block of any part of this document is permitted for non-profit educational purposes only. Copyright © 2015-2019 EDVOTEK, Inc., all rights reserved. 345-1905(1.1)

EDVOTEK

12

Preparing Students for PCR:

PCR is a ubiquitous technique that has become a cornerstone of modern biotechnology. Despite this, many of the details of the PCR process can be obscured behind the relatively simple protocol. Ensure that your students fully understand the science behind PCR by carefully reviewing the background information or incorporating supplemental lessons. Additional PCR lesson plans can be found at: www.edvotek.com/site/pdf/PCR_Lesson_Plan.pdf

Before You Start

- Picking up where you left off:
- If the samples were frozen at the end of Module I, they can be thawed on ice, or briefly held between the student's fingers.
 - Remind students that these are their extracted DNA samples which will be used for the PCR.
 - Review the steps of the protocol, ensure that each student or pair of students has the necessary reagents.
 - Ensure that the thermal cycler is programmed with the correct protocol and is functioning properly.

Pre-Lab Preparation

Detailed pre-lab preparation instructions can be found on page 21. Remember to plan for a stopping point if you will not be continuing directly to Module III.

IMPORTANT:
Prepare the controls before beginning the student protocol (see page 21). Controls should be run alongside student samples if possible.

Strategy & Planning

A properly mixed PCR reaction will be a light orange color after step 4. This is a good checkpoint to ensure that students have added the template and primers in the correct volumes.

Once the reaction mixture has been prepared (after step 4) we recommend proceeding through the PCR before stopping the experiment at the end of the module.


PCR is incredibly sensitive - each student must use a fresh pipet tip when mixing their samples.

Where's my DNA?


Experiment Overview



Where's my DNA?



Isolation of DNA from Cheek Cells




- Tip #1: Make sure the pellet is visible
- Tip #2: Resuspend VERY well in lysis buffer
- Tip #3: Water bath needs to be boiling

Isolation of DNA from Cheek Cells




- Tip #1: Run the Complete Control
- Tip #2: Students can use gloves to transfer PCR beads
- Tip #3: PCR product should be light orange

Isolation of DNA from Cheek Cells



- Tip #1: Perform digestion with the control
- Tip #2: Incubate for the whole 30 minutes
- Tip #3: Tubes can be incubated for longer than 30 min

Isolation of DNA from Cheek Cells




- Tip #1: Use Sybr Safe
- Tip #2: Use the DNA with the kit
- Tip #3: Load the Complete Control in lane 2
- Tip #4: Watch the gel closely - bands are small

Edvotek PLTW Experiments

Kit	Description
118-PLTW	Hypercholesterolemia Kit
225-PLTW	
235-PLTW	
268-PLTW	
301-PLTW	
303 - PLTW	
339 - PLTW	
990-PLTW	
116	
330	
345	
951	
953	Bioremediation by Oil-eating Bacteria Kit
956	

Mistakes ~~Hurt!~~ Happen!

Call



1-800-EDVOTEK

Online




YouTube
Facebook
Twitter
Pinterest

Email



info@Edvotek.com

Website




www.edvotek.com

We are ready to help!

Mistakes ~~Hurt!~~ Happen!

Call



1-800-EDVOTEK



Monday – Friday
8AM-5:30PM EST

We are ready to help!



What's Next?

For Orders & Technical Service:

- 1-800-EDVOTEK (1-800-338-6538)
- Web site: www.edvotek.com
- Email:
 - General - info@edvotek.com
 - Kelly - kellyb@edvotek.com

Check our YouTube Channel for instructional videos:

- www.youtube.com/edvotekinc

Like us on Facebook, Instagram, & Twitter for special offers & product updates: @edvotek

