



1121 5<sup>th</sup> Street NW  
Washington, DC 20001  
[www.edvotek.com](http://www.edvotek.com)

## Safety Data Sheet

Version 2  
Revision Date 08/19/2013

### 1. PRODUCT AND COMPANY IDENTIFICATION

Product name:	Calcium Chloride	Telephone:	1-800-338-6835
CaCl <sub>2</sub>			or 202-370-1500
Product Number:		Fax:	202-370-1501
Brand:	Edvotek	Emergency Phone #:	202-370-1500
Supplier:	Edvotek		
	1121 5 <sup>th</sup> Street NW		
	Washington, DC		
	20001		

### 2. HAZARDS IDENTIFICATION

#### Emergency Overview

#### OSHA Hazards:

Irritant.

#### GHS Classification

Acute toxicity (Category 5)

Eye Irritation (Category 2A)

#### Signal Word

Irritant

#### Pictogram



#### Hazard statement(s)

**H303** May be harmful if swallowed

Causes serious eye irritation

#### Precautionary statement(s)

**P305+P351+P338F** IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

#### HMIS Classification:

Health hazard: 2

Flammability: 0

Physical hazard: 0

**NFPA Rating:**

Health hazard: 2  
Fire: 0  
Reactivity Hazard: 0

**Potential Health Effects:**

**Inhalation** May be harmful if inhaled.

Causes respiratory tract irritation

**Skin** May be harmful if absorbed through skin. Causes skin irritation.

**Eyes** Causes eye irritation.

**Ingestion** May be harmful if swallowed.

### 3. COMPOSITION/INFORMATION ON INGREDIENTS

Formula: CaCl<sub>2</sub>

Molecular Weight: 110.98 g/mol

Chemical name	EC No.	CAS –No	Concentration
Calcium chloride	233-140-8	1004352-4	-

### 4. FIRST AID MEASURES

**General information**

No special measures required.

**If inhaled**

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult doctor in case of complaints.

**In case of skin contact**

Immediately wash off with soap and plenty of water. Get medical attention if symptoms develop.

**In case of eye contact**

Flush eyes with water as a precaution. Get medical attention immediately if irritation persists.

**If swallowed**

Never give anything by mouth to an unconscious person. Rinse mouth with water. Get medical attention if symptoms develop.

## 5. FIREFIGHTING MEASURES

### Conditions of flammability

Not flammable or combustible.

### Suitable extinguishing media

Use water spray, alcohol resistant foam, dry chemical or carbon dioxide.

### Special protective equipment for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

### Hazardous combustion products

Hazardous decomposition products formed under fire conditions. Hydrogen chloride gas, Calcium oxide

## 6. ACCIDENTAL RELEASE MEASURES

### Personal precautions

Prevent skin/eye contact. Use personal protective equipment as needed.

### Environmental precautions

Dilute with plenty of water. Do not allow to enter sewers/ surface or ground water.

### Methods and materials for containment and cleaning up

Clean up spills immediately. Use absorbent materials if needed and dispose into a chemical waste container.

## 7. HANDLING AND STORAGE

### Precautions for safe handling

No special measures necessary. Good laboratory technique should be used when handling.

### Conditions for safe storage

No special measures necessary. Store at room temperature.



Ignition temperature	no data available
Autoignition Temperature	no data available
Lower explosion limit	no data available
Upper explosion limit	no data available
Vapour pressure	no data available
Density (water=1)	no data available
Water solubility	no data available
Partition coefficient:	
n-octanol/water	no data available
Relative vapour Density	no data available
Odour Threshold	no data available
Evaporation rate	no data available

## 10. STABILITY AND REACTIVITY

### Chemical stability

Stable under recommended storage conditions.

oxide, Methyl vinyl ether, Calcium chloride is attacked by bromine Trifluoride

### Possibility of hazardous reactions

No further relevant data available.

### Hazardous decomposition products

Hazardous decomposition products formed under fire conditions.

### Conditions to avoid

No further relevant data available.

-Hydrogen chloride gas, Calcium oxide

**Materials to avoid** Strong acids, Borane/boron oxides, Zinc, Calcium

## 11. TOXICOLOGICAL INFORMATION

### Acute toxicity LD/LC50 values

no data available

### Respiratory or skin sensitization

no data available

### Serious eye damage/eye irritation

no data available

### Germ cell mutagenicity

no data available

**Skin corrosion/irritation**

no data available

**Carcinogenicity**

No components of this product are listed by IARC, ACGIH, NTP and OSHA.

**Reproductive toxicity**

no data available

**Aspiration hazard**

no data available

**Teratogenicity**

no data available

**Synergistic effects**

no data available

**12. ECOLOGICAL INFORMATION****Toxicity**

no data available

**Mobility in soil**

no data available

**Persistence and degradability**

no data available

**PBT and vPvB assessment**

no data available

**Bioaccumulative potential**

no data available

**Other adverse effects**

no data available

**13. DISPOSAL CONSIDERATIONS****Product**

Dispose of unused product, residues, and containers according to local, regional, state, and national regulations.

**Contaminated packaging**

Dispose of as unused product.

## 14. TRANSPORT INFORMATION

### **DOT (US)**

Not dangerous goods

### **IMDG**

Not dangerous goods

### **IATA**

Not dangerous goods

## 15. REGULATORY INFORMATION

### **OSHA Hazards**

Irritant

### **SARA 302 Components**

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

### **SARA 313 Components**

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the reporting levels established by SARA Title III, Section 313.

### **SARA 311/312 Hazards**

No SARA Hazards.

## **16. OTHER INFORMATION**

For educational use only. Not intended for human or animal diagnostic or therapeutic uses.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.