



1121 5th Street NW
 Washington, DC 20001
www.edvotek.com

Safety Data Sheet

Version 2
 Revision Date 10/29/2014

1. PRODUCT AND COMPANY IDENTIFICATION

| | | | |
|-----------------|--------------------------------|--------------------|-----------------|
| Product name: | Insect Cell Media | Telephone: | 1-800-338-6835 |
| Product Number: | | | or 202-370-1500 |
| Brand: | Edvotek | Fax: | 202-370-1501 |
| Supplier: | Edvotek | Emergency Phone #: | 202-370-1500 |
| | 1121 5 th Street NW | | |
| | Washington, DC | | |
| | 20001 | | |

2. HAZARDS IDENTIFICATION

Emergency Overview

Acute toxicity, harmful if swallowed.
 may cause an allergic skin reaction,
 may cause allergy or asthma
 symptoms or breathing difficulties if
 inhaled, suspected of damaging fertility
 or the unborn child.

GHS Classification

Respiratory sensitisation 1
 Skin sensitisation 1
 Reproductive toxicity 2

Signal Word

Danger

Pictogram



Hazard(s) Statements

H317 May cause an allergic skin
 reaction.
H334 May cause allergy or asthma
 symptoms
H361 Suspected of damaging fertility
 or the unborn child.

Precautionary statement(s)

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P280 Wear protective gloves.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.

P308 + P313 IF exposed or concerned: Get medical advice/ attention.

P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.

P363 Wash contaminated clothing before reuse.

P501 Dispose of contents/ container to an approved waste disposal plant

HMIS Classification:

Health hazard: 2

Flammability: 0

Physical hazard: 0

NFPA Rating:

Health hazard: 2

Fire: 0

Reactivity Hazard: 0

Potential Health Effects:

Inhalation May be harmful if inhaled.

Causes respiratory tract irritation

Skin May be harmful if absorbed through skin. Cause skin irritation.

Eyes May cause eye irritation.

Ingestion Harmful if swallowed.

4. COMPOSITION/INFORMATION ON INGREDIENTS

Mixture of the substances listed bellow with nonhazardous additions.

| Chemical name | EC No. | CAS –No | Concentration |
|-----------------------------|--------|---------|---------------|
| Insect Cell Media | | | 94% |
| Fetal Bovine Serum | | | 1% |
| Penicillin/ Streptomycin | | | 5% |

5. FIRST AID MEASURES

General information

No special measures required.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult doctor in case of complaints.

In case of skin contact

Immediately wash off with soap and plenty of water. Get medical attention if symptoms develop.

In case of eye contact

Flush eyes with water as a precaution. Get medical attention immediately if irritation persists.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Get medical attention if symptoms develop.

6. FIREFIGHTING MEASURES

Conditions of flammability

Not flammable or combustible.

Suitable extinguishing media

Use water spray, alcohol resistant foam, dry chemical or carbon dioxide.

Special protective equipment for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

Hazardous combustion products

Carbon oxides, nitrogen oxides (NO_x), Sulphur oxides

7. ACCIDENTAL RELEASE MEASURES

Personal precautions

Prevent skin/eye contact. Use personal protective equipment as needed.

Environmental precautions

Dilute with plenty of water. Do not allow to enter sewers/ surface or ground water.

Methods and materials for containment and cleaning up

Clean up spills immediately. Use absorbent materials if needed and dispose into a chemical waste container.

8. HANDLING AND STORAGE

Precautions for safe handling

No special measures necessary. Good laboratory technique should be used when handling.

Conditions for safe storage

No special measures necessary. Store at room temperature.

Incompatible Products

No known

8. EXPOSURE CONTROLS/ PERSONAL PROTECTION

Personal protective equipment

Respiratory protection

Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Wear protective gloves. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices.

Eye protection

Avoid eye contact with this material. Wear chemical safety goggles.

Skin and body protection

Use chemical gloves as needed. Use clean protective body covering clothing as needed to minimize contact with clothing and skin.

Hygiene measures

Handle in accordance with good industrial hygiene and safety practice. Wash and dry hands at the end of workday.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance

Form liquid
Colour colorless
Odor no data available

Safety data

| | |
|------------------------------|-------------------|
| pH | no data available |
| Melting point/freezing point | no data available |
| Boiling point | no data available |
| Flash point | no data available |
| Ignition temperature | no data available |
| Autoignition Temperature | no data available |
| Lower explosion limit | no data available |
| Upper explosion limit | no data available |
| Vapour pressure | no data available |
| Density (water=1) | no data available |
| Water solubility | no data available |
| Partition coefficient: | |
| n-octanol/water | no data available |
| Relative vapour Density | no data available |
| Odour Threshold | no data available |
| Evaporation rate | no data available |

10. STABILITY AND REACTIVITY

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

No further relevant data available.

Conditions to avoid

No further relevant data available.

Materials to avoid

No further relevant data available.

Hazardous decomposition products

No dangerous decomposition product

11. TOXICOLOGICAL INFORMATION

Acute toxicity LD/LC50 values

no data available

Serious eye damage/eye irritation

no data available

Respiratory or skin sensitization

no data available

Germ cell mutagenicity

no data available

Skin corrosion/irritation

no data available

Carcinogenicity

No components of this product are listed by IARC, ACGIH, NTP and OSHA.

Reproductive toxicity

no data available

Teratogenicity

no data available

Aspiration hazard

no data available

Synergistic effects

no data available

12. ECOLOGICAL INFORMATION**Toxicity**

no data available

Mobility in soil

no data available

Persistence and degradability

no data available

PBT and vPvB assessment

no data available

Bioaccumulative potential

no data available

Other adverse effects

no data available

13. DISPOSAL CONSIDERATIONS**Product**

Dispose of unused product, residues, and containers according to local, regional, state, and national regulations.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION

DOT (US)

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

15. REGULATORY INFORMATION

OSHA Hazards

Acute toxicity, harmful if swallowed, may cause an allergic skin reaction, may cause allergy or asthma symptoms or breathing difficulties if inhaled, suspected of damaging fertility or the unborn child.

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards.

16. OTHER INFORMATION

For educational use only. Not intended for human or animal diagnostic or therapeutic uses.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.