

Catalog #5019

Digital Orbital Shaker



Features:

- Digital speed control
- Speed range: 50-250 RPM
- Orbit 20 mm
- Digital timer: 1 min. to 99 h. 59 min.
- Maximum load capacity: 2.5 kg
- Orbital Shaker can be used in a cold room or incubator
- Temperature range: 5° to 40° C
- Power supply: 50 W

Our OR100 Digital Orbital Shaker is a powerful, variable-speed shaker which provides efficient orbital motion. Speed and time are under microprocessor control. The shaker can be applicable in microbiology, chemistry, immunology, biochemistry, and molecular biology laboratories.



1.800.EDVOTEK

| www.edvotek.com

| info@edvotek.com