



1121 5th Street NW
 Washington, DC 20001
www.edvotek.com

Safety Data Sheet

Version 2
 Revision Date 05/22/17

1. PRODUCT AND COMPANY IDENTIFICATION

| | | | |
|-----------------|--|--------------------|-----------------------------------|
| Product name: | TBE Electrophoresis Buffer Powder | Telephone: | 1-800-338-6835 or 202-370-1500 |
| Product Number: | | Fax: | 202-370-1501 |
| Brand: | Edvotek | Emergency Phone #: | 202-370-1500 |
| Supplier: | Edvotek 1121 5 th Street NW Washington, DC 20001 | | |

2. HAZARDS IDENTIFICATION

Emergency Overview

GHS Classification

Skin irritation (2)
 Eye irritation (2A)
 Reproductive toxicity (1B)
 Specific target organ toxicity - single exposure (3)
 Respiratory system, Acute aquatic toxicity (3)
 Chronic aquatic toxicity (3)

Signal Word

Danger

Pictogram



Hazard(s) Statements

H315 Causes skin irritation.
H319 Causes serious eye irritation.
H335 May cause respiratory irritation.
H360 . May damage fertility or the unborn child.
H412 Harmful to aquatic life with long lasting effects.

Precautionary statement(s)

P201 Obtain special instructions before use.

P202 Do not handle until all safety precautions have been read and understood.

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P264 Wash skin thoroughly after handling

P271 Use only outdoors or in a well-ventilated area.

P273 Avoid release to the environment

P280 Wear protective gloves/ eye protection/ face protection

P302 + P352 IF ON SKIN: Wash with plenty of soap and water

P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308 + P313: IF exposed or concerned: Get medical advice/ attention.

P321 Specific treatment (see supplemental first aid instructions on this label).

P332 + P313 If skin irritation occurs: Get medical advice/ attention.

P337 + P313 If eye irritation persists: Get medical advice/ attention

P362 Take off contaminated clothing and wash before reuse.

P403 + P233 Store in a well-ventilated place. Keep container tightly closed

P405 Store locked up.

P501 Dispose of contents/ container to an approved waste disposal plant

HMIS Classification:

Health hazard: 2

Flammability: 0

Physical hazard: 0

NFPA Rating:

Health hazard: 2

Fire: 0

Reactivity Hazard: 0

Potential Health Effects:

Inhalation May cause respiratory irritation.

Skin Cause skin irritation.

Eyes Causes serious eye irritation

Ingestion Harmful if swallowed.

3. COMPOSITION/INFORMATION ON INGREDIENTS

| Chemical name | EC No. | CAS –No | Concentration |
|---------------|-----------|------------|---------------|
| Boric Acid | 233-139-2 | 10043-35-3 | <50% |

4. FIRST AID MEASURES

General information

No special measures required.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult doctor in case of complaints.

In case of skin contact

Immediately wash off with soap and plenty of water. Get medical attention if symptoms develop.

In case of eye contact

Flush eyes with water as a precaution. Get medical attention immediately if irritation persists.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Get medical attention if symptoms develop.

5. FIREFIGHTING MEASURES

Conditions of flammability

Not flammable or combustible.

Suitable extinguishing media

Use water spray, alcohol resistant foam, dry chemical or carbon dioxide.

Special protective equipment for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

Hazardous combustion products

Carbon oxides, nitrogen oxides (NO_x), Borane/boron oxides, Sodium oxides

6. ACCIDENTAL RELEASE MEASURES

Personal precautions

Prevent skin/eye contact. Use personal protective equipment as needed.

Environmental precautions

Dilute with plenty of water. Do not allow to enter sewers/ surface or ground water.

Methods and materials for containment and cleaning up

Clean up spills immediately. Use absorbent materials if needed and dispose into a chemical waste container.

7. HANDLING AND STORAGE

Precautions for safe handling

Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Good laboratory technique should be used when handling.

Conditions for safe storage

No special measures necessary. Store at room temperature.

Incompatible Products

No known products

8. EXPOSURE CONTROLS/ PERSONAL PROTECTION

Personal protective equipment

Respiratory protection

Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Wear protective gloves. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices.

Eye protection

Avoid eye contact with this material. Wear chemical safety goggles.

Skin and body protection

Use chemical gloves as needed. Use clean protective body covering clothing as needed to minimize contact with clothing and skin.

Hygiene measures

Handle in accordance with good industrial hygiene and safety practice. Wash and dry hands at the end of workday.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance

Form powder
Colour white
Odor no data available

Safety data

| | |
|------------------------------|-------------------|
| pH | no data available |
| Melting point/freezing point | no data available |
| Boiling point | no data available |
| Flash point | no data available |
| Ignition temperature | no data available |
| Autoignition Temperature | no data available |
| Lower explosion limit | no data available |
| Upper explosion limit | no data available |
| Vapour pressure | no data available |
| Density (water=1) | no data available |
| Water solubility | no data available |
| Partition coefficient: | |
| n-octanol/water | no data available |
| Relative vapour Density | no data available |
| Odour Threshold | no data available |
| Evaporation rate | no data available |

10. STABILITY AND REACTIVITY

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

No further relevant data available.

Conditions to avoid

No further relevant data available.

Materials to avoid

No further relevant data available.

Hazardous decomposition products

No dangerous decomposition products

11. TOXICOLOGICAL INFORMATION

Acute toxicity LD/LC50 values

no data available

Serious eye damage/eye irritation

no data available

Respiratory or skin sensitization

no data available

Germ cell mutagenicity

no data available

Skin corrosion/irritation

no data available

Carcinogenicity

No components of this product are listed by IARC, ACGIH, NTP and OSHA.

Reproductive toxicity

no data available

Teratogenicity

no data available

Aspiration hazard

no data available

Synergistic effects

no data available

12. ECOLOGICAL INFORMATION**Toxicity**

no data available

Mobility in soil

no data available

Persistence and degradability

no data available

PBT and vPvB assessment

no data available

Bioaccumulative potential

no data available

Other adverse effects

no data available

13. DISPOSAL CONSIDERATIONS**Product**

Dispose of unused product, residues, and containers according to local, regional, state, and national regulations.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION**DOT (US)**

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

15. REGULATORY INFORMATION**OSHA Hazards**

Skin irritation

Eye irritation

Reproductive toxicity

Specific target organ toxicity - single exposure

Respiratory system, Acute aquatic toxicity

Chronic aquatic toxicity

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards.

16. OTHER INFORMATION

For educational use only. Not intended for human or animal diagnostic or therapeutic uses.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.