



1121 5th Street NW
Washington, DC 20001
www.edvotek.com

Safety Data Sheet

Version 2
Revision Date 09/19/2013

1. PRODUCT AND COMPANY IDENTIFICATION

Product name:	EDTA Solution	Telephone:	1-800-338-6835
Product Number:			or 202-370-1500
Brand:	Edvotek	Fax:	202-370-1501
Supplier:	Edvotek	Emergency Phone #:	202-370-1500
	1121 5 th Street NW		
	Washington, DC		
	20001		

2. HAZARDS IDENTIFICATION

Emergency Overview

GHS Classification

Eye irritation (Category 2A)

Signal Word

Irritant

Pictogram



Hazard(s) Statements

H319 Causes serious eye irritation

Precautionary statement(s)

P305 + P351 + P338 IF IN EYES:

Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

HMIS Classification:

Health hazard: 2

Flammability: 0

Physical hazard: 0

NFPA Rating:

Health hazard: 2
 Fire: 0
 Reactivity Hazard: 0

Potential Health Effects:

Inhalation May be harmful if inhaled.
 Causes respiratory tract irritation
Skin May be harmful if absorbed through skin. cause skin irritation.
Eyes May cause eye irritation.
Ingestion May be harmful if swallowed.

8. COMPOSITION/INFORMATION ON INGREDIENTS

Mixture of the substances listed bellow with nonhazardous additions.

Chemical name	EC No.	CAS –No	Concentration
EDTA	200-449-4	60-00-4	--
Water			

4. FIRST AID MEASURES

General information

No special measures required.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult doctor in case of complaints.

In case of skin contact

Immediately wash off with soap and plenty of water. Get medical attention if symptoms develop.

In case of eye contact

Flush eyes with water as a precaution. Get medical attention immediately if irritation persists.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Get medical attention if symptoms develop.

5. FIREFIGHTING MEASURES

Conditions of flammability

Not flammable or combustible.

Suitable extinguishing media

Use water spray, alcohol resistant foam, dry chemical or carbon dioxide.

Special protective equipment for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

Hazardous combustion products

Carbon oxides, Nitrogen oxides (NOx).

6. ACCIDENTAL RELEASE MEASURES

Personal precautions

Prevent skin/eye contact. Use personal protective equipment as needed.

Environmental precautions

Dilute with plenty of water. Do not allow to enter sewers/ surface or ground water.

Methods and materials for containment and cleaning up

Clean up spills immediately. Use absorbent materials if needed and dispose into a chemical waste container.

7. HANDLING AND STORAGE

Precautions for safe handling

No special measures necessary. Good laboratory technique should be used when handling.

Conditions for safe storage

No special measures necessary. Store at room temperature.

Incompatible Products

No known

Ignition temperature	no data available
Autoignition Temperature	no data available
Lower explosion limit	no data available
Upper explosion limit	no data available
Vapour pressure	no data available
Density (water=1)	no data available
Water solubility	no data available
Partition coefficient: n-octanol/water	no data available
Relative vapour Density	no data available
Odour Threshold	no data available
Evaporation rate	no data available

10. STABILITY AND REACTIVITY

Chemical stability

Stable under recommended storage conditions.

Conditions to avoid

No further relevant data available.

Possibility of hazardous reactions

No further relevant data available.

Materials to avoid

No further relevant data available.

Hazardous decomposition products

No dangerous decomposition products

11. TOXICOLOGICAL INFORMATION

Acute toxicity LD/LC50 values

no data available

Serious eye damage/eye irritation

no data available

Respiratory or skin sensitization

no data available

Germ cell mutagenicity

no data available

Skin corrosion/irritation

no data available

Carcinogenicity

No components of this product are listed by IARC, ACGIH, NTP and OSHA.

Reproductive toxicity

no data available

Teratogenicity

no data available

Aspiration hazard

no data available

Synergistic effects

no data available

12. ECOLOGICAL INFORMATION**Toxicity**

no data available

Mobility in soil

no data available

Persistence and degradability

no data available

PBT and vPvB assessment

no data available

Bioaccumulative potential

no data available

Other adverse effects

no data available

13. DISPOSAL CONSIDERATIONS

Product

Dispose of unused product, residues, and containers according to local, regional, state, and national regulations.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION

DOT (US)

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

15. REGULATORY INFORMATION

OSHA Hazards

No known OSHA hazards.

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards.

16. OTHER INFORMATION

For educational use only. Not intended for human or animal diagnostic or therapeutic uses.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.