



1121 5th Street NW
 Washington, DC 20001
www.edvotek.com

Safety Data Sheet

Version 2
 Revision Date 10/22/2015

1. PRODUCT AND COMPANY IDENTIFICATION

Product name: Wright stain	Telephone:	1-800-338-6835
Product Number:		or 202-370-1500
Brand: Edvotek	Fax:	202-370-1501
Supplier: Edvotek	Emergency Phone #:	202-370-1500
1121 5 th Street NW		
Washington, DC		
20001		

2. HAZARDS IDENTIFICATION

Emergency Overview

GHS Classification

Flammable liquids 2
 Acute toxicity, Oral 3
 Acute toxicity, Inhalation 3
 Acute toxicity, Dermal 3
 Eye irritation 2A
 Specific target organ toxicity - single exposure 1

Signal Word

Danger

Pictogram



Hazard statement(s)

H225 Highly flammable liquid and vapour.
H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled
H319 Causes serious eye irritation.
H370 Causes damage to organs.

Precautionary statement(s)

P210 Keep away from heat/sparks/open flames/hot surfaces.
- No smoking.

P233 Keep container tightly closed.

P240 Ground/bond container and receiving equipment.

P241 Use explosion-proof electrical/ventilating/ lighting/ equipment.

P242 Use only non-sparking tools.

P243 Take precautionary measures against static discharge.

P260 Do not breathe dust/ fume/ gas/ mist/ vapours/ spray.

P264 Wash skin thoroughly after handling.

P270 Do not eat, drink or smoke when using this product.

P271 Use only outdoors or in a well-ventilated area.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P310 + P330 IF

SWALLOWED: Immediately call a POISON CENTER or doctor/ physician. Rinse mouth.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

P304 + P340 + P311 IF INHALED:

Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER or doctor/ physician.

P305 + P351 + P338 IF IN EYES:

Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P307 + P311 IF exposed: Call a POISON CENTER or doctor/ physician.

P337 + P313 If eye irritation persists: Get medical advice/ attention.

P362 Take off contaminated clothing and wash before reuse.

P370 + P378 In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish.

P403 + P233 Store in a well-ventilated place. Keep container tightly closed.

P403 + P235 Store in a well-ventilated place.

HMIS Rating

Health hazard: 2

Chronic Health Hazard: 0

Flammability: 3

Physical Hazard 0

NFPA Rating

Health hazard: 2

Fire Hazard: 3

Reactivity Hazard: 0

3. COMPOSITION/INFORMATION ON INGREDIENTS

Synonyms: Eosin Methylene blue according to Wright

Chemical name	EC No.	CAS -No	Concentration
Methylthioninium chloride	200-515-2	61-73-4	>= 10 - < 20 %
Disodium 2-(2,4,5,7-tetrabromo-6-oxido-3-oxoxanthen-9-yl)benzoate	241-409-6	17372-87-1	>= 10 - < 20 %
Methanol	200-659-6	67-56-1	>= 70 - < 90 %

4. FIRST AID MEASURES

General information

No special measures required.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult doctor in case of complaints.

In case of skin contact

Immediately wash off with soap and plenty of water. Get medical attention if symptoms develop.

In case of eye contact

Flush eyes with water as a precaution. Get medical attention immediately if irritation persists.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Get medical attention if symptoms develop.

5. FIREFIGHTING MEASURES

Conditions of flammability

Not flammable or combustible.

Suitable extinguishing media

Use water spray, alcohol resistant foam, dry chemical or carbon dioxide.

Special protective equipment for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

Hazardous combustion products

Carbon oxides, nitrogen oxides (NO_x), Sulphur oxides, Hydrogen chloride gas, Hydrogen bromide gas, Sodium Oxides

6. ACCIDENTAL RELEASE MEASURES

Personal precautions

Prevent skin/eye contact. Use personal protective equipment as needed.

Environmental precautions

Dilute with plenty of water. Do not allow to enter sewers/ surface or ground water.

Methods and materials for containment and cleaning up

Clean up spills immediately. Use absorbent materials if needed and dispose into a chemical waste container.

7. HANDLING AND STORAGE

Precautions for safe handling

Avoid contact with skin and eyes.
Avoid inhalation of vapour or mist.
Use explosion-proof equipment. Keep away from sources of ignition - No smoking.

Conditions for safe storage

Keep container tightly closed in a dry and well-ventilated place. Containers which are opened must be carefully resealed and kept upright to prevent leakage.

Incompatible Products

No known products

8. EXPOSURE CONTROLS/ PERSONAL PROTECTION

Personal protective equipment

Respiratory protection

Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Wear protective gloves. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices

Eye protection

Avoid eye contact with this material. Wear chemical safety goggles.

Skin and body protection

Use chemical gloves as needed. Use clean protective body covering clothing as needed to minimize contact with clothing and skin.

Hygiene measures

Handle in accordance with good industrial hygiene and safety practice. Wash and dry hands at the end of workday.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance

Form	liquid
Color	blue
Odor	no data available

Safety data

pH	no data available
Melting point/freezing point	no data available
Boiling point	no data available
Flash point	no data available
Ignition temperature	no data available
Autoignition Temperature	no data available
Lower explosion limit	no data available
Upper explosion limit	no data available
Vapor pressure	no data available
Density (water=1)	no data available
Water solubility	no data available
Partition coefficient: n-octanol/water	no data available
Relative vapour Density	no data available
Odour Threshold	no data available
Evaporation rate	no data available

10. STABILITY AND REACTIVITY

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

No further relevant data available

Conditions to avoid

No data available

Materials to avoid

No further relevant data available

Hazardous decomposition products

No known products

11. TOXICOLOGICAL INFORMATION**Acute toxicity LD/LC50 values**

no data available

Serious eye damage/eye irritation

no data available

Respiratory or skin sensitization

no data available

Germ cell mutagenicity

no data available

Skin corrosion/irritation

no data available

Carcinogenicity

No components of this product are listed by IARC, ACGIH, NTP and OSHA.

Reproductive toxicity

no data available

Teratogenicity

no data available

Aspiration hazard

no data available

Synergistic effects

no data available

12. ECOLOGICAL INFORMATION**Toxicity**

no data available

Mobility in soil

no data available

Persistence and degradability

no data available

PBT and vPvB assessment

no data available

Bioaccumulative potential

no data available

Other adverse effects

no data available

13. DISPOSAL CONSIDERATIONS**Product**

Dispose of unused product, residues, and containers according to local, regional, state, and national regulations.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION**DOT (US)**

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

15. REGULATORY INFORMATION

OSHA Hazards

Flammable liquids 2

Acute toxicity, Oral 3

Acute toxicity, Inhalation 3

Acute toxicity, Dermal 3

Eye irritation 2A

Specific target organ toxicity - single exposure 1

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards.

16. OTHER INFORMATION

For educational use only. Not intended for human or animal diagnostic or therapeutic uses.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.