

CATALOG #546

EDVOTEK® Incubation Oven



Specifications:

Temperature Range (H): 5°C above ambient to 60°C

Temperature Accuracy: 0.5°C (at 37°C)

Temperature Uniformity: +/- 1.5°C (at 37°C)

Temperature Increment: 1°C

Dimensions (w x d x h):

Exterior: 13.2 x 14.5 x 18.7 in. / 33.5 x 37 x 47.5 cm

Interior: 10.3 x 9.3 x 12.8 in. / 26 x 23.5 x 32.5 cm

Capacity: 20L (0.75 ft³)

Weight: 15 lbs / 6.5 kg

Warranty: 2 Years

Electrical: 100 to 240V, 50-60Hz

With digital temperature control, the EDVOTEK® Incubation Oven provides convenient “set and walk away” operation, eliminating the need for external thermometers and repetitive “fine tuning” of an analog control knob. Simply choose the desired temperature and the incubator immediately begins to heat up while accurately monitoring the chamber temperature and conveniently displaying the temperature in real time on the large LED control panel.



Despite a modest footprint of 13.2 x 14.5", our incubator features a large internal chamber, capable of storing flasks and bottles up to 2 liters. In addition, our incubator includes two adjustable shelves and an internal outlet for powering small shakers, rockers and other instrumentation.



EDVOTEK®

1.800.EDVOTEK
www.edvotek.com