

Catalog #552

White LED Transilluminator

Features:

- 25 x 25 cm Viewing Area
- 90 Lumens LED Chips
- Power Rating: 8 Watts
- Universal Voltage: 110/220 V
- Frequency: 50/60 Hz
- Two-Year Warranty
- Made in USA



Our White LED Transilluminator features a spacious 25 x 25 cm viewing area illuminated by long life LEDs and is housed in a slim aluminum body. It's designed to safely enhance the visualization of DNA stained with FlashBlue™, proteins stained with Coomassie Blue and autoradiograms.

CARE AND MAINTENANCE

- Do not expose the transilluminator to organic solvents or immerse it in water.
- Clean the viewing surface with a damp sponge or paper towel.

1.800.EDVOTEK
www.edvotek.com

EDVOTEK®