

Luria Broth Media (100 grams)

Catalog #611

EXCITE | EXPLORE | ENGAGE



Storage: Room Temperature

The media contains hydrolyzed casein, yeast extract and salts for the growth of bacteria. Twenty-five grams is added to every liter of distilled water. The liquid media is then autoclaved at 121° C for 15 minutes.

EDVOTEK products are intended for educational use only. They are not to be used for diagnostic or drug purposes, nor administered to or consumed by humans or animals.

Conditions: EDVOTEK, Inc. warrants the materials it sells to meet EDVOTEK's performance specifications in effect at the time of shipment. EDVOTEK, Inc.'s liability under this warranty is limited to the repair of defects or the replacement of the product at its option. It is subject to receipt of reasonable proof by the customer that the defect is embraced within the terms of the warranty. All claims made under warranty must be presented to EDVOTEK, Inc. within 30 days following the date of delivery of the product to customer. EDVOTEK, Inc. does not accept responsibility for user conversion, abusive use, lack of, or improper maintenance of equipment.

1.800.EDVOTEK | www.edvotek.com | info@edvotek.com