


# X-Gal

## Catalog #614

EXCITE | EXPLORE | ENGAGE



**Storage: Below 0°C**

Synonym: 5-Bromo-4-chloro-3-indolyl  $\beta$ -D-galactopyranoside  
Formula Weight: 408.6

Conditions: EDVOTEK, Inc. warrants the materials it sells to meet EDVOTEK's performance specifications in effect at the time of shipment. EDVOTEK, Inc.'s liability under this warranty is limited to the repair of defects or the replacement of the product at its option. It is subject to receipt of reasonable proof by the customer that the defect is embraced within the terms of the warranty. All claims made under warranty must be presented to EDVOTEK, Inc. within 30 days following the date of delivery of the product to customer. EDVOTEK, Inc. does not accept responsibility for user conversion, abusive use, lack of, or improper maintenance of equipment.

1.800.EDVOTEK | [www.edvotek.com](http://www.edvotek.com) | [info@edvotek.com](mailto:info@edvotek.com)