



1121 5th Street NW
 Washington, DC 20001
www.edvotek.com

Safety Data Sheet

Version 3
 Revision Date 05/14/2025

PRODUCT AND COMPANY IDENTIFICATION

Product name: Nickel Solution	Telephone: 1-800-338-6835
Product Number:	or 202-370-1500
Brand: Edvotek	Fax: 202-370-1501
Supplier: Edvotek	Emergency Phone #: 202-370-1500
1121 5 th Street NW	
Washington, DC 20001	

HAZARDS IDENTIFICATION

Emergency Overview

OSHA Hazards:

Acute toxicity Oral, Acute toxicity Inhalation, ^[L]_[SEP] Skin irritation Respiratory sensitization, Skin sensitization, ^[L]_[SEP] Germ cell mutagenicity, ^[L]_[SEP] Carcinogenicity, Reproductive toxicity ^[L]_[SEP] Specific target organ toxicity , Inhalation Acute aquatic toxicity, Chronic aquatic toxicity.

Skin irritation 2
 Respiratory sensitization 1
 Skin sensitization 1
 Germ cell mutagenicity 2
 Carcinogenicity 1A
 Reproductive toxicity 1B^[L]_[SEP]
 Specific target organ toxicity 1
 Inhalation Acute aquatic toxicity 1
 Chronic aquatic toxicity 1

GHS Classification

Acute toxicity Oral 4
 Acute toxicity Inhalation 4

Signal Word

Danger

Pictogram



Hazard statement(s)

H302 + H332 Harmful if swallowed or if inhaled

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

H341 Suspected of causing genetic defects.

H350 May cause cancer

H360 May damage fertility or the unborn child.

H372 Causes damage to organs through prolonged or repeated exposure if inhaled

H410 Very toxic to aquatic life with long lasting effects.

Precautionary statement(s)

P201 Obtain special instructions before use

P202 Do not handle until all safety precautions have been read and understood.

P264 Wash skin thoroughly after handling

P270 Do not eat, drink or smoke when using this product.

P272 Contaminated work clothing should not be allowed out of the workplace

P273 Avoid release to the environment

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P312 IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. Rinse mouth.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P308 + P313 IF exposed or concerned: Get medical advice/ attention.

P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.

P362 Take off contaminated clothing and wash before reuse.

P391 Collect spillage.

P501 Dispose of contents/ container to an approved waste disposal plant.

HMIS Classification:

Health hazard: 2

Flammability: 0

Physical hazard: 0

NFPA Rating:

Health hazard: 2

Fire: 0

Reactivity Hazard: 0

Potential Health Effects:

Inhalation: May be harmful if inhaled.

Skin: May cause skin irritation.

Eyes: May cause eye irritation.

Ingestion: Harmful if swallowed.

COMPOSITION/INFORMATION ON INGREDIENTS

Formula: NiO4S · 6H2O

Molecular weight: 262.85 g/mol

Chemical name	EC No.	CAS -No	Concentration
Nickel sulphate hexahydrate	232-104-9 _{SEP}	10101-97-0 _{SEP}	15 mg/ml

FIRST AID MEASURES

General information

No special measures required.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult doctor in case of complaints.

In case of skin contact

Immediately wash off with soap and plenty of water. Get medical attention if symptoms develop.

In case of eye contact

Flush eyes with water as a precaution. Get medical attention immediately if irritation persists.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Get medical attention if symptoms develop.

FIREFIGHTING MEASURES

Conditions of flammability

Not flammable or combustible.

Suitable extinguishing media

Use water spray, alcohol resistant foam, dry chemical or carbon dioxide.

Special protective equipment for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

Hazardous combustion products

No known products

ACCIDENTAL RELEASE MEASURES**Personal precautions**

Prevent skin/eye contact. Use personal protective equipment as needed.

Environmental precautions

Dilute with plenty of water. Do not allow to enter sewers/ surface or ground water.

Methods and materials for containment and cleaning up

Clean up spills immediately. Use absorbent materials if needed and dispose into a chemical waste container.

HANDLING AND STORAGE**Precautions for safe handling**

No special measures necessary. Good laboratory technique should be used when handling.

Conditions for safe storage

No special measures necessary. Store at room temperature.

Incompatible Products

No known products

8. EXPOSURE CONTROLS/ PERSONAL PROTECTION

Personal protective equipment

Respiratory protection

Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Wear protective gloves. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices.

Eye protection

Avoid eye contact with this material. Wear chemical safety goggles.

Skin and body protection

Use chemical gloves as needed. Use clean protective body covering clothing as needed to minimize contact with clothing and skin.

Hygiene measures

Handle in accordance with good industrial hygiene and safety practice. Wash and dry hands at the end of workday.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance

Form	liquid
Color	green
Odor	no data available

Safety data

pH	no data available
Melting point/freezing point	no data available

Boiling point	no data available
Flash point	no data available
Ignition temperature	no data available
Autoignition Temperature	no data available
Lower explosion limit	no data available
Upper explosion limit	no data available
Vapor pressure	no data available
Density (water=1)	no data available
Water solubility	no data available
Partition coefficient:	
n-octanol/water	no data available
Relative vapour Density	no data available
Odour Threshold	no data available
Evaporation rate	no data available

10. STABILITY AND REACTIVITY

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

No further relevant data available.

Conditions to avoid

No further relevant data available.

Materials to avoid

No further relevant data available.

Hazardous decomposition products

No dangerous decomposition products known.

TOXICOLOGICAL INFORMATION

Acute toxicity LD/LC50 values

no data available

Serious eye damage/eye irritation

no data available

Respiratory or skin sensitization

no data available

Germ cell mutagenicity

no data available

Skin corrosion/irritation

no data available

Carcinogenicity

No components of this product are listed by IARC, ACGIH, NTP and OSHA.

Reproductive toxicity

no data available

Teratogenicity

no data available

Aspiration hazard

no data available

Synergistic effects

no data available

12. ECOLOGICAL INFORMATION**Toxicity**

no data available

Mobility in soil

no data available

Persistence and degradability

no data available

PBT and vPvB assessment

no data available

Bioaccumulative potential

no data available

Other adverse effects

no data available

13. DISPOSAL CONSIDERATIONS

Product

Dispose of unused product, residues, and containers according to local, regional, state, and national regulations.

Contaminated packaging

14. TRANSPORT INFORMATION

Dispose of as unused product.

DOT (US)

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

15. REGULATORY INFORMATION

OSHA Hazards

Acute toxicity Oral, Acute toxicity Inhalation, ^[]_[SEP] Skin irritation Respiratory sensitization, Skin sensitization, ^[]_[SEP] Germ cell mutagenicity, ^[]_[SEP] Carcinogenicity, Reproductive toxicity ^[]_[SEP] Specific target organ toxicity , Inhalation Acute aquatic toxicity, Chronic aquatic toxicity.

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards.

16. OTHER INFORMATION

For educational use only. Not intended for human or animal diagnostic or therapeutic uses.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.